



# 7710

## Bed scale approved

**SOEHNLE**  
PROFESSIONAL

- S** Easy-lift system simplifies the lifting of bed wheels
- S** Protection class IP65 for use in damp locations

- > Approved to class III, Medical device directive (MPG) class 1
- > Easy-lift system: for easy lifting
- > Suitable for bed wheels with approx. 130 - 160 mm diameter
- > Distance piece for smaller bed wheels can be ordered optionally
- > Alphanumeric keyboard for the input of organizational data (e.g. bed weights)
- > RS 232 interface for data transfer to a computer or printer
- > BMI, hold and dialysis function
- > Plus-minus weighing
- > Fixed tare value memory (up to 100 bed weights), manual input of tare values
- > LCD digit height: 20 mm
- > Units in kg or lb
- > Protection class IP65



Order number	7710.01.001	7710.02.001
Version	with trolley	without trolley
Maximum capacity	500 kg	500 kg
Division	100 g / 200 g	100 g / 200 g
Dimensions	875 x 470 x 1,080 mm	
Scale weight	33 kg	25 kg

### Accessories

Data cable 2 m (display - PC/IT)

Order number 2550.03.002



## Bed scales

for hospitals and  
rehabilitation clinics



# 7711

## Bed scale Bluetooth™ approved

**SOEHNE**  
**PROFESSIONAL**  
**PROFESSIONAL**

- S** Cordless data transfer via Bluetooth™
- S** Easy-lift system simplifies the lifting of bed wheels

- › Approved to class III, Medical device directive (MPG) class 1
- › Cordless data transfer via Bluetooth™ (max. operating distance 10 m)
- › Easy-lift system: for easy lifting
- › Suitable for bed wheels with approx. 130 - 160 mm diameter
- › Distance piece for smaller bed wheels can be ordered optionally
- › Alphanumeric keyboard for the input of organizational data (e.g. bed weights)
- › BMI, hold and dialysis function
- › Plus-minus weighing
- › Fixed tare value memory (up to 100 bed weights), manual input of tare values
- › RS 232 interface for data transfer to a computer or printer
- › Units in kg or lb
- › Display with mains operation
- › Foot sensors with rechargeable battery operation (continuous operating time 10 h)
- › 2 battery packs, charging and storage station (included in scope of delivery)



Order number	7711.01.001	7711.02.001
Version	with trolley	without trolley
Maximum capacity	500 kg	500 kg
Division	100 g / 200 g	100 g / 200 g
Dimensions	875x470x1,080 mm	
Scale weight	40 kg	22.8 kg

### Accessories

Data cable 2 m (display - PC/IT)  
Order number 2550.03.002

Retrofit kit for wheel diameter 100 - 120 mm  
Order number 2503.04.001

Retrofit kit for wheel diameter 90 - 100 mm  
Order number 2503.04.002



Medguard  
Professional Healthcare Supplies  
Tel: +353 (01) 835 2411 Email: service@medguard.ie  
Fax: +353 (01) 835 2461 Web: www.medguard.ie