

# 7801

## Personal scale with case

- S** *Ideal for mobile use*
- S** *Weighing platform with separate display*

- › Fast 1-button operation
- › Easy to read LCD display, 20 mm digit height
- › Tare function for mother-child-weighing
- › Hold function for better readability
- › Model without carry case, includes wall bracket
- › Display made of high quality plastic
- › Flat weighing platform with easy-to-clean surface
- › Cable length of approx. 2 m from display to platform
- › Units in kg or lb
- › Energy-saving auto-off function
- › Rechargeable battery operation with integrated charger (operating time 50 h, charging time 5.5 h)



### Without carry case

Order number	7801.75.002
Maximum capacity	200 kg
Division	100 g
Dimensions	Weighing platform: 310x275x34 mm
Scale weight	7 kg

### With carry case

Order number	7801.70.002
Maximum capacity	200 kg
Division	100 g
Dimensions	Weighing platform: 310x275x34 mm
Scale weight	9 kg

Medguard  
Professional Healthcare Supplies  
Tel: +353 (01) 835 2411 Email: [service@medguard.ie](mailto:service@medguard.ie)  
Fax: +353 (01) 835 2461 Web: [www.medguard.ie](http://www.medguard.ie)

