


Emergency Trolley (CI-09114, CI-09116)

- Made mainly of ABS Plastic, Aluminum and Stainless Steel
- Light and easy to clean
- Weight: 36 kg approx.



Medguard Professional Healthcare

Tel: +353 (0) 1 835 2411

Fax: +353 (0) 1 969 5009

Email: service@medguard.ie

Website: www.medguard.ie

Specifications	Emergency Trolley (CI-09114, CI-09116)
Dimensions	620 x 480 x (H) 970 mm
Drawer Size (inside) (3)	428 x 338 x (H) 50 mm
Drawer Size (inside) (2)	425 x 335 x (H) 105 mm
Transparent Cardiac Board	Yes
Built-in Perimetral Protection	Yes
Castors	Ø 75 mm
Side Waste Container and Hanging Box	Yes
Rotating Defibrillator Holder Tray	Yes
Oxygen Tank Holder	Yes
IV Pole	Yes
Side Straight Handle	Yes
Sliding Handwriting Board	Yes
Stainless Steel Guard Rail	Yes
Power Outlets and Line (240V)	Yes

Medguard Professional Healthcare

Tel: +353 (0) 1 835 2411

Fax: +353 (0) 1 969 5009

Email: service@medguard.ie

Website: www.medguard.ie