

Professional Healthcare Supplies

PROBE'S **EB02.002**









It is a large isothermal ampoule holder that allows the user to carry a high number of ampoules and medicines. Thanks to its 4 detachable trays with identifying tags, the health professional will reduce the time of location and administration of the medicines, thus achieving a better patient assistance.

*1 reusable cold gel included.

Features

- Night/day reflective pullers (allows a quicker opening of the bag).
- Padded handle.
- Cardholder on the front of the bag.
- Mesh side pockets & Isothermal material.
- Probe's has two detachable trays to hold a total of 68 small units, 52 medium sizes, 368 big units and 2 extra large. It has also an internal pocket of extra material and 4 mesh pockets to carry some cold gels.
- Removable middle compartment with transparent window for syringes, rubber tourniquet, gauzes, etc.
- 4 mesh pockets on the sides to place the cold gels.

Technical data

- Measures: 28 x 21 x 9 cm
- Capacity: 5,29 L
- **Weight**: 0.65 kg
- Material: 600 D Polyester. Washable
- Colour: Red
- Maximum recommended contents weight: 3,5 kg

Medguard Professional Healthcare