

ERKA.

THE ORIGINAL

Made in Germany since 1889



Medguard Professional Healthcare

Tel: +353 (0) 1 835 2411

Email: service@medguard.ie

Fax: +353 (0) 1 969 5009

Website: www.medguard.ie

Product catalogue 2013



Welcome to the family.

Welcome to **ERKA**. We are a traditional family company that has been developing outstanding medical equipment of the highest precision and quality with great passion for more than 120 years. Among professional users, our name stands for groundbreaking innovation, continuity and reliability. As the pioneer of mechanical blood pressure measurement, our devices have proven themselves millions of times around the world and set the standard for all future developments.

With kind regards,

B. Kallweitz



Content

ERKA. The Original. Perfection made in Germany. Since 1889.	2
Our Service	8
Contact	9
Mein ERKA.	10
Products	11
Sphygmomanometers Professional	11
Perfect Aneroid	
Profi	
Kobold	
Vario	
E3000	
Klinik ERKA.meter	
Sphygmomanometers ERKA. SWITCH 2.0	25
SWITCH 2.0 Komfort	
SWITCH 2.0 Simplex	
Sphygmomanometers Homecare	32
ERKA.test	
Simplitest	
Kobitest	
Stethoscopes Dual Membrane	38
Precise	
Sensitive	
Classic	
Stethoscopes ERKA.Membrane	43
Finesse ²	
Finesse	
ERKA.phon	
ERKA.phon duo	
High Pressure Cuffs	50
SWITCH 2.0 HPC	
Cuffs	53
Green Cuff	
Classic Cuff	
Spare parts	59

Your satisfaction is our mission.

Our claim at **ERKA** is to provide medical professionals with the most precise and efficient equipment possible, meeting their needs to pursue their profession. Our development efforts focus on ease of operation, durability and functionality. We do not view innovation as technology confusion but as a mechanical masterpiece. Which is one more reason that every product with the **ERKA** logo also bears the "Made in Germany" quality seal.



Origins are the future.

ERKA. makes history and you can help write it in the future: The company was founded in Berlin by Richard Kallmeyer in the year 1889. With the development of the first industrially produced blood pressure measurement device, the first milestone in the company history followed in 1914. In 1927, the "Erkameter" became the global epitome of blood pressure measurement. While **ERKA.** distributes its products in more than 100 countries today, production remains in just one country – Germany – as it has been from day one.



You cannot win trust.

You have to earn it. Through continuity and sustainability. As an owner-managed family company in the fourth generation, we at **ERKA** place great emphasis on long-term customer relationships defined by mutual esteem. We view our customers as permanent members of our family, and feel a great commitment to meeting their professional challenges. While this has a lot to do with responsibility, we are motivated first and foremost by enthusiasm for the cause. Your cause.



Performance thrives on passion.

The special quality of **ERKA** products arises from the identification of all our employees with their tasks and their clients, dedicating all of their experience and passion to the individual needs of the customer. Long-term employment relationships with low turnover play a role along with a high level of professionalism and the service orientation of every individual. Before and after the sale. At **ERKA**, we are always here for you, no matter where you are.



Made in Germany, made for you.

In Germany the saying goes "Qualität beginnt mit Qual" – quality begins with an ordeal. What does this mean for **ERKA**? Before a "Perfect Aneroid" blood pressure measurement device for example is delivered, it has to pass a total of 45 individual product inspections in our Quality Management Check (QMC) – from the material inspection to the detailed inspection to the comprehensive verification of functionality. We are aware of our responsibility and refuse to rely on others, ensuring that chance cannot impair the quality of your work and results in the future.



Efficiency demands excellence.

For **ERKA** products, the best is only just good enough. From the quality of the base materials to the selection of suppliers, from friendly service and product usability to the environmental compatibility of packaging and environmentally friendly production processes, we highly value excellence in all respects. While we certainly do not offer the cheapest products in the market, they are ultimately the most advantageous for our customers. Efficiency pays off over the long term.



ERKA 360° Service. Everything in one place.

What **ERKA**. product is right for me? How do I use my **ERKA**. product correctly? What laws and directives do I have to observe when operating my **ERKA**. product? Who helps me with calibration and medical-technical inspections? Where can I find an **ERKA**. service partner near me?

At **ERKA**. we offer full service for our products. Either at service@erka.org or in cooperation with one of our many partners!



Talk to us. We are happy to help you.

Thank you!



Focus on you.

Take advantage of the opportunity to add your **ERKA** medical product to our database and benefit from our many services.

We take the initiative to inform you of calibration due dates and make recommendations for the timely replacement of wear parts. You couldn't ask for an easier or more efficient partnership!



Professionalism according to ERKA.:

All **ERKA**. blood pressure measurement equipment complies with the international EN 1060 standard and Medical Device Directive 93/42/EEC (MDD). Quality is guaranteed by our quality management system certified according to DIN EN ISO. FDA approval has also been obtained for select models.

ERKA. blood pressure measurement devices are impact-resistant based on recognised standards such as ANSI/AAMI and SP 10, and meet the requirements for use in an MR room.

Especially hardened copper-beryllium membranes ensure a long lifespan with consistently high precision. Clinically tested, latex-free and with an accuracy of 1% across the entire measurement range, **ERKA**. sets the standard in blood pressure measurement.

Sphygmomanometers | Professional



ERKA Perfect Aneroid

The precision and durability standard.

Proven in practice millions of times worldwide, ERKA has been setting the standard in blood pressure measurement for decades with the Perfect-Aneroid.

This unique success is based on the ERKA Quality Management Check, a series of 45 inspections during production, including a computer-controlled 100% final inspection. This allows us to guarantee that the user receives the best possible blood pressure measurement devices for their profession, and nothing less – this is the claim supported by our daily work!



The scales optimised for legibility are available in 48 and 56 mm diameters. Right and left-handed models are offered.



The clinic model with pressure gauge protection cover for challenging conditions in everyday clinical use.



Maximum operator comfort thanks to the continuously variable automatic valve.



High-gloss, polished chromium-plated housing for ultimate surface hygiene, ergonomically optimised design, ultra-flat construction and sensitive threaded discharge valve – the perfection of blood pressure measurement.



The selection of ERKA Green Cuff products is comprehensive. They are durable, skin-friendly, biocompatible and stand for ultimate wearer comfort.

	Cuff · Type · Color	Size 1: 10-15 cm	Size 2: 14-21 cm	Size 3: 20,5-28 cm	Size 4: 27-35 cm	Size 5: 34-43 cm	Size 6: 42-54 cm
Type: Standard Scale: 48 mm	superb · Rapid · green	201.28482	201.27482	201.26482	201.20482	201.24482	201.29482
	superb · Hook · green	-	-	201.26282	201.20282	201.24282	-
	superb · D-Ring · green	-	-	201.46882	201.40882	201.44882	-
Type: Standard Scale: 56 mm	superb · Rapid · green	203.28482	203.27482	203.26482	203.20482	203.24482	203.29482
	superb · Hook · green	-	-	203.26282	203.20282	203.24282	-
	superb · D-Ring · green	-	-	203.46882	203.40882	203.44882	-
Type: Automatic Scale: 48 mm	superb · Rapid · green	202.28482	202.27482	202.26482	202.20482	202.24482	202.29482
	superb · Hook · green	-	-	202.26282	202.20282	202.24282	-
	superb · D-Ring · green	-	-	202.46882	202.40882	202.44882	-
Type: Automatic Scale: 56 mm	superb · Rapid · green	204.28482	204.27482	204.26482	204.20482	204.24482	204.29482
	superb · Hook · green	-	-	204.26282	204.20282	204.24282	-
	superb · D-Ring · green	-	-	204.46882	204.40882	204.44882	-
Type: Clinic Scale: 48 mm	superb · Rapid · green	215.28482	215.27482	215.26482	215.20482	215.24482	215.29482
	superb · Hook · green	-	-	215.26282	215.20282	215.24282	-
	superb · D-Ring · green	-	-	215.46882	215.40882	215.44882	-
Type: Clinic Scale: 56 mm	superb · Rapid · green	214.28482	214.27482	214.26482	214.20482	214.24482	214.29482
	superb · Hook · green	-	-	214.26282	214.20282	214.24282	-
	superb · D-Ring · green	-	-	214.46882	214.40882	214.44882	-

Left handed models are available without extra charge. Just note on your order.

	Cuff · Type · Color	Size 1: 10-15 cm	Size 2: 14-21 cm	Size 3: 20,5-28 cm	Size 4: 27-35 cm	Size 5: 34-43 cm	Size 6: 42-54 cm
Type: Standard Scale: 48 mm	superb · Rapid · green	271.28482	271.27482	271.26482	271.27482	271.24482	271.29482
	superb · Hook · green	-	-	271.26282	271.27282	271.24282	-
	superb · D-Ring · green	-	-	271.46882	271.40882	271.44882	-
Type: Standard Scale: 56 mm	superb · Rapid · green	273.28482	273.27482	273.26482	273.27482	273.24482	273.29482
	superb · Hook · green	-	-	273.26282	273.27282	273.24282	-
	superb · D-Ring · green	-	-	273.46882	273.40882	273.44882	-
Type: Automatic Scale: 48 mm	superb · Rapid · green	272.28482	272.27482	272.26482	272.27482	272.24482	272.29482
	superb · Hook · green	-	-	272.26282	272.27282	272.24282	-
	superb · D-Ring · green	-	-	272.46882	272.40882	272.44882	-
Type: Automatic Scale: 56 mm	superb · Rapid · green	274.28482	274.27482	274.26482	274.27482	274.24482	274.29482
	superb · Hook · green	-	-	274.26282	274.27282	274.24282	-
	superb · D-Ring · green	-	-	274.46882	274.40882	274.44882	-
Type: Clinic Scale: 48 mm	superb · Rapid · green	275.28482	275.27482	275.26482	275.27482	275.24482	275.29482
	superb · Hook · green	-	-	275.26282	275.27282	275.24282	-
	superb · D-Ring · green	-	-	275.46882	275.40882	275.44882	-
Type: Clinic Scale: 56 mm	superb · Rapid · green	276.28482	276.27482	276.26482	276.27482	276.24482	276.29482
	superb · Hook · green	-	-	276.26282	276.27282	276.24282	-
	superb · D-Ring · green	-	-	276.46882	276.40882	276.44882	-

Left handed models are available without extra charge. Just note on your order.

ERKA. Profi

Precision by professionals, for professionals.

Conceived by **ERKA.** with a focus on the cost-benefit aspect, this blood pressure measurement device in a class of its own combines the efficient use of high-quality materials with maximum measurement accuracy and ultimate ease of use. For example, the **ERKA.** Profi uses dual-hose membrane protection technology for lasting precision, while numerous cuffs that can be combined permit versatile application in everyday hospital use. Of course the **ERKA.** Quality Management Check with its numerous material inspections and functional checks guarantees precision and exact measurements for years with this instrument as well!



The scales optimised for legibility are available in 48 and 56 mm diameters.



The clinic model with pressure gauge protection cover for challenging conditions in everyday clinical use.



ERKA. Profi – impact-resistant based on internationally recognised standards such as ANSI/AAMI SP 10.



Ergonomic and balanced, made of lightweight materials and equipped with a sensitive threaded discharge valve, the **ERKA.** Profi guarantees maximum ease of use.



The selection of **ERKA.** Green Cuff products is comprehensive. They are durable, skin-friendly, biocompatible and stand for ultimate wearer comfort.

	Cuff · Type · Color	Size 1: 10-15 cm	Size 2: 14-21 cm	Size 3: 20,5-28 cm	Size 4: 27-35 cm	Size 5: 34-43 cm	Size 6: 42-54 cm
Type: Standard Scale: 48 mm	superb · Rapid · grey	228.28492	228.27492	228.26492	228.20492	228.24492	228.29492
	superb · Hook · grey	-	-	228.26292	228.20292	228.24292	-
	superb · D-Ring · grey	-	-	228.46892	228.40892	228.44892	-
	smart · Rapid · grey	228.38392	228.37392	228.36392	228.30392	228.34392	228.39392
Type: Standard Scale: 56 mm	superb · Rapid · grey	229.28492	229.27492	229.26492	229.20492	229.24492	229.29492
	superb · Hook · grey	-	-	229.26292	229.20292	229.24292	-
	superb · D-Ring · grey	-	-	229.46892	229.40892	229.44892	-
	smart · Rapid · grey	229.38392	229.37392	229.36392	229.30392	229.34392	229.39392
Type: Klinik Scale: 48 mm	superb · Rapid · grey	230.28492	230.27492	230.26492	230.20492	230.24492	230.29492
	superb · Hook · grey	-	-	230.26292	230.20292	230.24292	-
	superb · D-Ring · grey	-	-	230.46892	230.40892	230.44892	-
	smart · Rapid · grey	230.38392	230.37392	230.36392	230.30392	230.34392	230.39392

ERKA. Kobold

Reliable and exact in everyday medical application.

The ERKA. Kobold is a robust and reliable partner for all professional users.

Thanks to the ergonomic, balanced design combined with optimisation for low weight and the highly sensitive threaded discharge valve, maximum ease of use is guaranteed in the daily routine of the medical professional – for every measurement.



Anodised aluminium housing in five different colours for optimised surface hygiene and cleaning results.



The ERKA. Kobold. Available with three different cuff sets.



Anodised aluminium housing in five different colours for optimised surface hygiene and cleaning results.



The robust scale with a diameter of 48 mm is optimised for legibility.



The selection of ERKA. Green Cuff products is comprehensive. They are durable, skin-friendly, biocompatible and stand for ultimate wearer comfort.

Cuff · Type · Color		Size 1: 10-15 cm	Size 2: 14-21 cm	Size 3: 20,5-28 cm	Size 4: 27-35 cm	Size 5: 34-43 cm	Size 6: 42-54 cm
Type: Standard							
Color:							
Black: Colorcode 223	superb · Rapid · grey	Colorcode.28493	Colorcode.27493	Colorcode.26493	Colorcode.20493	Colorcode.24493	Colorcode.29493
Red: Colorcode 224	superb · Hook · grey	-	-	Colorcode.26293	Colorcode.20293	Colorcode.24293	-
Blue: Colorcode 225	superb · D-Ring · grey	-	-	Colorcode.46893	Colorcode.40893	Colorcode.44893	-
Gold: Colorcode 226	smart · Rapid · grey	Colorcode.38393	Colorcode.37393	Colorcode.36393	Colorcode.30393	Colorcode.34393	Colorcode.39393
Silver: Colorcode 227							
Scale: 48 mm							
Type: Standard							
Scale: 48 mm							
with Cuff-Set 1:	superb · Rapid · grey	227.71393	Set component	Set component	-	-	-
includes sizes 1-3	smart · Rapid · grey	227.74393	Set component	Set component	-	-	-
Color: Silver							
Type: Standard							
Scale: 48 mm							
with Cuff-Set 2:	superb · Rapid · grey	-	-	227.71293	Set component	Set component	-
includes sizes 3-5	superb · Hook · grey	-	-	227.72293	Set component	Set component	-
Color: Silver	superb · D-Ring · grey	-	-	227.74293	Set component	Set component	-
	smart · Rapid · grey	-	-	227.73293	Set component	Set component	-
Type: Standard							
Scale: 48 mm							
with Cuff-Set 3:	superb · Rapid · grey	227.71193	Set component	Set component	Set component	Set component	Set component
includes sizes 1-6	smart · Rapid · grey	227.73193	Set component	Set component	Set component	Set component	Set component
Color: Silver							

The standard set is equipped with a silver colored pressure gauge. Other colors are available on request.

ERKA. Vario

Diligence, comfort and flexibility.

The **ERKA. Vario** – a stationary blood pressure measurement device meeting all requirements of the modern medical professional.

The precision pressure gauge is recessed in a premium, impact-resistant, easy to clean ABS housing. The large cuff basket, the easy to use spiral hose, the inflation bulb with the precision adjustable threaded discharge valve and optimised readability on the oversize scale through anti-glare glass guarantee ultimate functionality.



ERKA. Vario base model with wall bracket, available separately.



ERKA. Vario base model with pedestal, available separately.



ERKA. Vario base model with rail holder, available separately.



Three different cuff sets, quick and simple to exchange with the **ERKA. quick-closure**, are available for optimum measuring.



The selection of **ERKA. Green Cuff** products is comprehensive. They are durable, skin-friendly, biocompatible and stand for ultimate wearer comfort.

	Cuff · Type · Color	Size 1: 10-15 cm	Size 2: 14-21 cm	Size 3: 20,5-28 cm	Size 4: 27-35 cm	Size 5: 34-43 cm	Size 6: 42-54 cm
ERKA Vario desktop model with one cuff	superb · Rapid · green	285.28481	285.27481	285.26481	285.20481	285.24481	285.29481
	superb · Hook · green	-	-	285.26281	285.20281	285.24281	-
	superb · D-Ring · green	-	-	285.46881	285.40881	285.44881	-
eco type:	smart · Rapid · grey	285.38391	285.37391	285.36391	285.30391	285.34391	285.39391
ERKA Vario desktop model with cuff-set 1: includes size 1-3	superb · Rapid · green	285.71381	Set component	Set component	-	-	-
	eco type:	smart · Rapid · grey	285.73391	Set component	Set component	-	-
ERKA Vario desktop model with cuff-set 2: includes size 3-5	superb · Rapid · green	-	-	285.71281	Set component	Set component	-
	superb · Hook · green	-	-	285.72281	Set component	Set component	-
	superb · D-Ring · green	-	-	285.74281	Set component	Set component	-
eco type:	smart · Rapid · grey	-	-	285.73291	Set component	Set component	-
ERKA Vario desktop model with cuff-set 3: includes size 1-6	superb · Rapid · grün	285.71181	Set component	Set component	Set component	Set component	Set component
	eco type:	smart · Rapid · grey	285.73191	Set component	Set component	Set component	Set component
ERKA. Vario, standard model	superb · Rapid · green	288.28481	288.27481	288.26481	288.20481	288.24481	288.29481
	superb · Hook · green	-	-	288.26281	288.20281	288.24281	-
	superb · D-Ring · green	-	-	288.46881	288.40881	288.44881	-
eco type:	superb · Rapid · grey	288.28491	288.27491	288.26491	288.20491	288.24491	288.29491
	smart · Rapid · grey	288.38391	288.37391	288.36391	288.30391	288.34391	288.39391

Equipment for the **ERKA.** Vario standard model


Stand	Enables the use as a stand, with 5 rolls	995.00003	All-purpose holder (height 15-50mm)	Enables the use as a rail model, sleuable	295.28700
Wall holder	Enables the use as a wall holder, sleuable	295.00003			

ERKA.meter 3000

Thoroughness in perfection.

The name **ERKA.meter** has represented the STANDARD and is the ORIGINAL in blood pressure measurement since 1927.

Available in two equipment versions, the **ERKA.meter 3000** is the latest generation of our most successful product family, combining a modern design with highly innovative technical details. Millions of medical professionals worldwide trust the **ERKA.meter**. The precision of an **ERKA.** quicksilver blood pressure measurement device remains unrivalled to this day.



Conceived as a rectangle, the shatterproof capillary tube guarantees that measurements are read precisely without parallax errors while simultaneously reducing quicksilver use by 40%, helping to protect the environment.



When the device is closed, the automatic protection feature fully drains the quicksilver from the capillary tube back into the reservoir and prevents return flow.



The robust ABS housing is sturdy, safe and impact-resistant.



The newly developed pressure relief valve prevents pressure in excess of 315 mmHg. Damage to the air filter leading to leaking quicksilver is excluded.



The selection of **ERKA.** Green Cuff products is comprehensive. They are durable, skin-friendly, biocompatible and stand for ultimate wearer comfort.

	Cuff · Type · Color	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6
Type: Standard	superb · Rapid · green	177.28481	177.27481	177.26481	177.20481	177.24481	177.29481
	superb · Hook · green	-	-	177.26281	177.20281	177.24281	-
	superb · D-Ring · green	-	-	177.46881	177.40881	177.44881	-
	smart · Rapid · green	177.38381	177.37381	177.36381	177.30381	177.34381	177.39381
Typ: eco	superb · Rapid · grey	176.28491	176.27491	176.26491	176.20491	176.24491	176.29491
	superb · Hook · grey	-	-	176.26291	176.20291	176.24291	-
	superb · D-Ring · grey	-	-	176.46891	176.40891	176.44891	-
	smart · Rapid · grey	176.38391	176.37391	176.36391	176.30391	176.34391	176.39391

Klinik **ERKA**.meter

A tradition of diligence.

The universal classic in the **ERKA**.meter family!

Blood pressure measurement with quicksilver sets the standard for precision and reliability. However, we as the manufacturer must guarantee the responsible use of quicksilver. Which is why the highest standards for quality and functionality apply to all blood pressure measurement equipment in the **ERKA**.meter family and in the course of our **ERKA**. Quality Management Check – for each and every device!



Klinik **ERKA**.meter base model with wall bracket, available separately.



Klinik **ERKA**.meter base model with pedestal, available separately.



Klinik **ERKA**.meter base model with rail holder, available separately.



Three different cuff sets, quick and simple to exchange with the **ERKA**. quick-closure, are available for optimum measuring.



The selection of **ERKA**. Green Cuff products is comprehensive. They are durable, skin-friendly, biocompatible and stand for ultimate wearer comfort.

	Cuff · Type · Color	Size 1: 10-15 cm	Size	Size 3: 20,5-28 cm	Size 4: 27-35 cm	Size 5: 34-43 cm	Size 6: 42-54 cm
ERKA .meter desk model with one cuff	superb · Rapid · green	109.28481	109.27481	109.26481	109.20481	109.24481	109.29481
	superb · Hook · green	-	-	109.26281	109.20281	109.24281	-
	superb · D-Ring · green	-	-	109.46881	109.40881	109.44881	-
eco type:	smart · Rapid · grey	109.38391	109.37391	109.36391	109.30391	109.34391	109.39391
ERKA .meter desk model with cuff-set 1: includes size 1-3	superb · Rapid · green	109.71381	Set component	Set component	-	-	-
	eco type:	smart · Rapid · grey	109.73391	Set component	Set component	-	-
ERKA .meter desk model with cuff-set 2: includes size 3-5	superb · Rapid · green	-	-	109.71281	Set component	Set component	-
	superb · Hook · green	-	-	109.72281	Set component	Set component	-
	superb · D-Ring · green	-	-	109.74281	Set component	Set component	-
eco type:	smart · Rapid · grey	-	-	109.73281	Set component	Set component	-
ERKA .meter desk model with cuff-set 3: includes size 1-6	superb · Rapid · green	109.71181	Set component	Set component	Set component	Set component	Set component
	eco type:	smart · Rapid · grey	109.73191	Set component	Set component	Set component	Set component
ERKA . Klink ERKA meter, standard model with one cuff	superb · Rapid · green	116.28481	116.27481	116.26481	116.20481	116.24481	116.29481
	superb · Hook · green	-	-	116.26281	116.20281	116.24281	-
	superb · D-Ring · green	-	-	116.46881	116.40881	116.44881	-
eco type:	smart · Rapid · grey	116.38391	116.37391	116.36391	116.30391	116.34391	116.39391

Equipment for the Klinik **ERKA**.meter standard model

(independently from the selected cuff)

Stand	Enables the use as a stand, with 5 rolls	995.00003	All-purpose holder (height 15-50mm)	Enables the use as a rail model, slewable	295.28700
Wall holder	Enables the use as a wall holder, slewable	295.00003			

ERKA. SWITCH 2.0

Maximum innovation based on 123 years of experience.

As the blood pressure measurement pioneer, ERKA. has developed a new generation of blood pressure measurement equipment inspired by practical experience and in close cooperation with users around the world.



Sphygmomanometers | ERKA SWITCH

The claim.

Economic efficiency and ecologically responsible action.


The implementation.

On the new **ERKA. SWITCH 2.0**, the load cell can be replaced intuitively by the user within seconds.


The result.

- Saves time for the medical technician since the old load cell is simply replaced by a new, pre-calibrated load cell as part of the metrology inspection
- Saves time for the user, since the device is once again ready for use on the patient within seconds
- Saves time for administration, since elaborate purchasing, logistics and warehousing processes are greatly simplified
- Conserves resources – when purchasing a replacement or performing the metrology inspection, only the load cell has to be replaced

ERKA Switch 2.0 Comfort



ERKA Switch 2.0 Simplex und HPC:



ERKA. SWITCH 2.0

The performance.

- Absolute accuracy guaranteed by the ERKA. Quality Management Check (QMC), including a computer-supported 100% final inspection and more
- Sophisticated devices with convincing everyday clinic performance due to impact resistance and MR compatibility
- Innovative solutions such as left-right switching for maximum ease of use



ERKA. SWITCH 2.0

The possibilities.

Thanks to its innovative design, the handling and functionality of the ERKA. SWITCH 2.0 complies with all international standards, making it the first global blood pressure measurement device – Made in Germany.

ERKA SWITCH 2.0



ERKA SWITCH 2.0 Comfort



ERKA SWITCH 2.0 Simplex



ERKA SWITCH 2.0 HPC

ERKA. Switch 2.0 Comfort

Efficiency down to the smallest detail.

ERKA. as the pioneer of blood pressure measurement developed the ERKA. SWITCH 2.0 Comfort – based on the ERKA. requirements for a modern blood pressure measurement device, and in close cooperation with medical professionals worldwide.

Outstanding in quality and precision, comfortable in function and handling, convincing in efficiency and benefits, responsible in the use of ecological resources – the modern standard defined by ERKA.



The requirements for use in an MR room are met.



The polycarbonate material is extremely easy to clean as well as impact-resistant based on internationally recognised standards such as ANSI/AAMI SP10. The scale with a diameter of 56 mm is optimised for legibility.



The selection of ERKA. Green Cuff products is comprehensive. They are durable, skin-friendly, biocompatible and stand for ultimate wearer comfort.



ERKA. Switch 2.0 Simplex

Elementary and precise.

Precise measurement and efficient handling – the essentials made by ERKA!

The **ERKA. SWITCH 2.0 Simplex** appeals with absolute precision guaranteed by the **ERKA. Quality Management Check** and the consistent realisation of our claim:

Minimum effort for daily handling with maximum ecological responsibility.



The requirements for use in an MR room are met.



Replacing the load cell? Child's play. Exchange it for a new, calibrated load cell in a few simple steps. Without the help of qualified medical-technical personnel!



The polycarbonate material is extremely easy to clean as well as impact-resistant based on internationally recognised standards such as ANSI/AAMI SP10. The scale with a diameter of 56 mm is optimised for legibility.



Ergonomic and balanced, made of lightweight materials and equipped with a sensitive threaded discharge valve, the **ERKA. SWITCH 2.0** guarantees maximum ease of use.

	Cuff · Type · Color	Size 1: 10-15 cm	Size 2: 14-21 cm	Size 3: 20,5-28 cm	Size 4: 27-35 cm	Size 5: 34-43 cm	Size 6: 42-54 cm
SWITCH 2.0 Simplex, one cuff	superb · D-Ring · grey	-	-	294.46891	294.40891	294.44891	-
SWITCH 2.0 Comfort, one cuff	superb · Rapid · grey	293.28493	293.27493	293.26493	293.20493	293.24493	293.29493
	superb · Hook · grey	-	-	293.26293	293.20293	293.24293	-
	superb · D-Ring · grey	-	-	293.46893	293.40893	293.44893	-
	smart · Rapid · grey	293.38393	293.37393	293.36393	293.30393	293.34393	293.39393
SWITCH 2.0 Comfort with cuff-set 1: Size 1-3	superb · Rapid · grey	293.71393	Part of the set	Part of the set	-	-	-
	smart · Rapid · grey	293.74393	Part of the set	Part of the set	-	-	-
SWITCH 2.0 Comfort with cuff-set 2: Size 3-5	superb · Rapid · grey	-	-	293.71293	Part of the set	Part of the set	-
	superb · Hook · grey	-	-	293.72293	Part of the set	Part of the set	-
	superb · D-Ring · grey	-	-	293.74293	Part of the set	Part of the set	-
	smart · Rapid · grey	-	-	293.73293	Part of the set	Part of the set	-
SWITCH 2.0 Comfort with cuff-set 3: Size 1-6	superb · Rapid · grey	293.71193	Part of the set	Part of the set	Part of the set	Part of the set	Part of the set
	smart · Rapid · grey	293.73193	Part of the set	Part of the set	Part of the set	Part of the set	Part of the set
	Size 1: 500 ml	Size 2: 1000 ml					
SWITCH 2.0 HPC Pressure infusion cuff	822.00000	832.00000					

Accessories for SWITCH 2.0

	Housing for Comfort	Housing for HPC and Simplex (with hook)	Hook (separat)	safety clip (separat)	Body with ball for Comfort
SWITCH 2.0	000.30200	000.30201	995.00139	995.00140	000.30100

ERKA. Homecare.

**Your doctor has been trusting
us since 1889 – and so should you!**

Millions of doctors worldwide have been putting their trust in **ERKA.** products for more than 120 years. They are convinced by the precision and quality of our measuring instruments. Every **ERKA.** blood pressure management device, whether for professionals or private users, passes a computer-supported 100% final inspection.

In cooperation with medical professionals around the world, we developed the familiar **ERKA.** blood pressure measurement devices for continuous blood pressure monitoring at home. If your



Sphygmomanometers | Homecare

ERKA.test

The perfection of personal monitoring and reliability.

Based on the **ERKA**. Perfect-Aneroid, the blood pressure measurement standard since decades!

Dual-hose membrane protection technology guarantees consistent precision, the impact-resistant design maximises suitability for everyday use and the high-gloss, polished, chromium-plated housing is free of pores and very easy to clean.



The scales optimised for legibility are available in 48 and 56 mm diameters.



The clinic model with pressure gauge protection cover for challenging conditions in everyday clinical use.



Maximum operator comfort thanks to the continuously variable automatic valve, available as an option.



High-gloss, polished chromium-plated housing for ultimate surface hygiene, ergonomically optimised design, ultra-flat construction and sensitive threaded discharge valve – the perfection of blood pressure measurement.



Ultra-flat, high-quality stethoscope conceived especially for blood pressure measurement.

ERKA. Simplitest

The homecare professional.

The ERKA. Simplitest is a premium blood pressure measurement device, combining the efficient use of high-quality materials with maximum measurement accuracy and ultimate ease of use.

Of course the ERKA. Quality Management Check with its numerous material inspections and functional checks guarantees precision and exact measurements for years with the ERKA. Simplitest as well!



The scales optimised for legibility are available in 48 and 56 mm diameters.



The clinic model with pressure gauge protection cover for challenging conditions in everyday use.



Impact-resistant pressure gauge based on internationally recognised standards such as ANSI/AAMI SP10.



Precise results by using a stethoscope conceived especially for blood pressure measurement.



The selection of ERKA. Green Cuff products is comprehensive. They are durable, skin-friendly, biocompatible and stand for ultimate wearer comfort.

ERKA. Kobitest

Reliable and robust.

The ERKA. Kobitest is a robust and reliable everyday partner.

Thanks to the ergonomic, balanced design combined with optimisation for low weight and the highly sensitive threaded discharge valve, maximum ease of use is guaranteed in the daily routine of the medical professional – for every measurement!



The scale fabricated with a diameter of 48 mm is optimised for legibility.



Precise results by using a stethoscope conceived especially for blood pressure measurement.



The ultra-flat stethoscope conceived especially for blood pressure measurement ensures precise measurement results.

	Cuff · Type · Color	Size 3	Size 4	Size 5
Type: ERKA.test Scale: 48 mm	superb · D-Ring-Homecare · green	205.06982	205.00982	205.04982
Type: ERKAtest Scale: 56 mm	superb · D-Ring-Homecare · green	206.06982	206.00982	206.04982
Type: ERKAtest, automatic Scale: 48 mm	superb · D-Ring-Homecare · green	216.06982	216.00982	216.04982
Type: ERKAtest, automatic Scale: 56 mm	superb · D-Ring-Homecare · green	217.06982	217.00982	217.04982
Type: Simplitest Scale: 48 mm	superb · D-Ring-Homecare · grey	240.06992	240.00992	240.04992
Type: Simplitest Scale: 56 mm	superb · D-Ring-Homecare · grey	218.06992	218.00992	218.04992
Type : Kobitest (Farbe: Schwarz) Scale : 48 mm	superb · D-Ring-Homecare · grey	219.06993	219.00993	219.04993
Type: Switch Homecare Scale: 56 mm	superb · D-Ring-Homecare · grey	292.06993	292.00993	292.04993

Replacement stethoscop:

	Material · Surface	black	grey	blue	green	burgundy	red	bright green	violet	dark grey
Type: flat	Brass · bright chrome	543.00010	543.00000	543.00020	543.00055	543.00060	543.00030	543.00051	543.00080	543.00005
Size: adult	Brass · matt chrome	543.01010	543.01000	543.01020	543.01055	543.01060	543.01030	543.01051	543.01080	543.01005
	Aluminium · anodized	544.00010	544.00000	544.00020	544.00055	544.00060	544.00030	544.00051	544.00080	544.00005

Further spare parts

Blood pressure journal w. case	Cuff D-Ring-Homecare		Size 3	Size 4	Size 5
000.04000.4270	Single-channel-tube cuff for Kobitest	green	025.46883	025.40883	025.44883
		grey	025.46893	025.40893	025.44893
	Double-channel tube cuff for ERKAtest, Simplitest	green	025.46881	025.40881	025.44881
		grey	025.46891	025.40891	025.44891

ERKA. Stethoskope. Excellent hearing.

We make all **ERKA.** stethoscopes of the highest quality brass – naturally made in Germany.

In addition to the high specific weight, our stethoscopes take advantage of another outstanding property of brass: in an elaborate process based on many years of experience accumulated by our employees, every **ERKA.** stethoscope receives a high-gloss polish and chromium plating free of pores. Compared to all other materials/surfaces, this process guarantees the best cleaning results and therefore also the best surface hygiene!



Stethoscopes

ERKA. Dual Membrane. Perfect hearing.

The optimum transmission and amplification of body sounds requires the thinnest possible membranes. Yet they also have to withstand the demands of everyday medical use.

With the patented ERKA. Dual Membrane, ERKA. offers a solution that effectively combines both requirements: the synthetic membrane facing the skin is resistant to daily stresses, protecting the highly sensitive and extra-thin copper-beryllium membrane underneath it. The space created by joining the two membranes acts as an additional resonance chamber and pre-amplifier – for perfect hearing.




ERKA. Precise

Hear everything.

ERKA. has redefined the term "high-end stethoscope" with the ERKA. Precise!

The open funnel side equipped with a ring as protection against cold allows the medical professional to hear lower frequencies in particular. The other side is equipped with the ERKA. Dual Membrane which, thanks to its unique, patented structure, permits the most sensitive auscultation of both low and high frequencies.



Composition and effect of the patented ERKA. Dual Membrane.



The floating membrane prevents sound transmission at deformation. A ring as protection against cold improves patient comfort and protects the membrane.



Two parallel interior channels run inside the preformed tube – available in numerous different colours. This avoids friction noise which occurs when two separate tubes are used.



Thanks to the spring integrated in the tube, acoustic signals are conducted from the bell to the ear with virtually no loss.




The ergonomically shaped, super soft ear buds offer a good fit for enhanced wearer comfort and prevent acoustic loss. Naturally, the soft ear buds shown are also included as an alternative.

ERKA. Sensitive

The auscultation standard.

The combination of the sound chamber geometry and ERKA. dual membrane technology permits optimum auscultation of both low and high tones. Second to none for blood pressure measurement and cardiology. The perfection of inner values is consummately complemented by the exterior design of this high-end stethoscope. High-gloss polished with chrome plating that is free of pores, the brass bell offers perfect surface hygiene.



Composition and effect of the patented ERKA. Dual Membrane:



The floating membrane prevents sound transmission at deformation. A ring as protection against cold improves patient comfort and protects the membrane.



Two parallel interior channels run inside the preformed tube – available in numerous different colours. This avoids friction noise which occurs when two separate tubes are used.



Thanks to the spring integrated in the tube, acoustic signals are conducted from the bell to the ear with virtually no loss.




The ergonomically shaped, super soft ear buds offer a good fit for enhanced wearer comfort and prevent acoustic loss. Naturally, the soft ear buds shown are also included as an alternative.

ERKA. Classic

Perceive the essentials.

The best of the ERKA. tradition is reduced to the essentials!

A sound chamber geometry that has been perfected over the course of decades in combination with the modern ERKA. Dual Membrane provides the best possible auscultation across the entire frequency range. Auscultation is optimised by the high specific weight of the brass body, reliably eliminating ambient noise.



Composition and effect of the patented ERKA. Dual Membrane:er patentierten ERKA. Dual Membrane.



The floating membrane prevents sound transmission at deformation. A ring as protection against cold improves patient comfort and protects the membrane.



The brass bell. Free of pores and chrome-plated in high-gloss polished or matt version.



Thanks to the spring integrated in the tube, acoustic signals are conducted from the bell to the ear with virtually no loss.





The ergonomically shaped, super soft ear buds offer a good fit for enhanced wearer comfort and prevent acoustic loss. Naturally, the soft ear buds shown are also included as an alternative.

Product overview

ERKA. Stethoscope with **ERKA.** Dual Membrane

	Material · Surface	black	blue	burgundy	grey	green	red
Sensitive	Brass · bright chrome	525.00000	525.00020	525.00060	525.00005	525.00055	525.00030
Precise	Brass · bright chrome	531.00000	531.00020	531.00060	531.00005	531.00055	531.00030
	Brass · bright chrome	531.01000	531.01020	531.01060	531.01005	531.01055	531.01030
Classic	Brass · bright chrome	570.00000	570.00020	570.00060	570.00005	570.00055	570.00030
	Brass · matt chrome	570.01000	570.01020	570.01060	570.01005	570.01055	570.01030

Special colors as pink  or bright blue  are available by request.

Spare-part overview

	Type	Sensitive	Precise	Classic	Finesse	Finesse ²	ERKA.phon duo	ERKA.phon
Membrane	Dual-Membrane	000.40200	000.40200	000.40200	-	-	-	-
	Standard Membrane infant	-	-	-	-	000.40102	-	-
	Standard Membrane chid	-	-	-	000.40100	000.40100	000.40100	000.40100
	Standard Membrane adult	-	-	-	000.40000	000.40000	000.40000	000.40000
non-chill ring	Dual-Membrane	000.53104	000.53104	000.53104	-	-	-	-
	Size infant	-	-	-	-	000.53103	-	-
	Size child	-	-	-	000.53142	000.53142	000.53142 [*]	000.53142 [*]
	Size adult	-	-	-	000.53141	000.53141	000.53141 [*]	000.53141 [*]
Non-chill ring for funnel		-	000.53304	-	000.53304	-	000.53304	-
Eartips	Soft	-	-	000.45004	-	-	000.45004	000.45004
	Soft	000.45104	000.45104	-	000.45104	000.45104	-	-
	Super Soft	000.44900	000.44900	000.44900	000.44900	000.44900	000.44900 [*]	000.44900 [*]

* optional for
ERKA.phon Stethoscops

ERKA. Membrane. **Universal hearing.**

The **ERKA.** Membrane is suitable for the auscultation of both high and low frequencies, making it the ideal and universal instrument for medical professionals.

The **ERKA.** Membrane is installed with convex basic stress, which is consistently maintained by the metal membrane holder ring. This innovative design noticeably improves auscultation results.



ERKA. Finesse

The professional standard.

The stethoscope for every medical professional!

With a switchable double bell, this stethoscope covers every eventuality of daily life in a medical practice. The open funnel side is intended especially for low frequencies. The membrane side covers the full frequency range. A chrome-plated, high-gloss polished surface is free of pores, and all of the materials used can be cleaned perfectly and quickly for maximum surface hygiene.



ERKA. standard membrane with convex tension for outstanding resonance and sound transmission. g für eine hervorragende Resonanz- und Schallübertragung.



The brass bell, polished or with a matt chrome finish – both versions available in adult and child sizes. A ring for protection against cold on both sides maximises patient comfort.



Two parallel interior channels run inside the preformed tube – available in numerous different colours. This avoids friction noise which occurs when two separate tubes are used.



Thanks to the spring integrated in the tube, acoustic signals are conducted from the bell to the ear with virtually no loss.



The ergonomically shaped, super soft ear buds offer a good fit for enhanced wearer comfort and prevent acoustic loss. Naturally, the soft ear buds shown are also included as an alternative.

ERKA. Finesse²

One for all.

Since every patient is unique, **ERKA.** makes a stethoscope that is helpful to the medical professional in virtually all patient situations.

The switchable double bell features two membrane sides, each of which makes the entire frequency range perfectly audible. Since the **ERKA.** standard membranes are installed with convex basic stress, they guarantee the best hearing. A ring as protection against cold is installed on each side to improve patient comfort.



The bells produced in three different sizes are combined with each other as follows: ADULT & CHILD as well as CHILD & INFANT



The brass bell is available in the following versions: free of pores and chrome-plated with a high-gloss polish, or chrome-plated and matt.



Two parallel interior channels run inside the preformed tube – available in numerous different colours. This avoids friction noise which occurs when two separate tubes are used.



The spring separately integrated in the tube ensures that acoustic signals are conducted to the ear with virtually no loss.



The ergonomically shaped, soft ear buds offer a good fit for enhanced wearer comfort and prevent acoustic loss. Naturally, the hard ear buds shown are also included as an alternative.

	Material · Surface	black	grey	blue	green	burgundy	red
Finesse	Brass · bright chrome	550.00000	550.00005	550.00020	550.00055	550.00060	550.00030
	Brass · matt chrome	550.01000	550.01005	550.01020	550.01055	550.01060	550.01030
Finesse Child	Brass · bright chrome	549.00000	549.00005	549.00020	549.00055	549.00060	549.00030
	Brass · matt chrome	549.01000	549.01005	549.01020	549.01055	549.01060	549.01030
Finesse ²	Brass · bright chrome	535.00000	535.00005	535.00020	535.00055	535.00060	535.00030
	Brass · matt chrome	535.01000	535.01005	535.01020	535.01055	535.01060	535.01030
Finesse ² Child	Brass · bright chrome	536.00000	536.00053	536.00020	536.00055	536.00060	536.00030
	Brass · matt chrome	536.01000	536.01005	536.01020	536.01055	536.01060	536.01030

Special colors as pink  or bright blue  are available by request.

Spare-part overview

[illegible]

ERKA.phon Duo

Straightforward. Universal.

The ERKA.phon product family offers an extensive selection of models for the various areas of application of a stethoscope.

Available for adults and children, the ERKA.phon Duo features a switchable double bell. The open funnel side is designed especially for low frequencies, the membrane side for the entire frequency range. The interplay of high-end components in the ERKA.phon series ensures the best possible acoustics.



Since the ERKA. standard membranes are installed with convex basic stress, they guarantee the best hearing. A ring as protection against cold is available for the membrane side as an option.



The bell – available in brass/ chrome-plated brass with high-gloss finish/ chrome-plated with matt finish or anodised aluminium – is quick and easy to clean.



The tube is available in nine colour variations.



To ensure that signals captured by the bell are transmitted to the ears without interference, the spring was not integrated in the tube duct.



The ERKA.phon Duo is also available in a special version for teaching applications.



ERKA.phon

Hearing blood pressure.

The ERKA.phon is the optimum stethoscope for blood pressure measurement, since the flat bell is easy to position underneath the cuff. In the interplay with the remaining high-quality components, ERKA. standard membranes installed with convex basic stress guarantee the best possible hearing across the entire frequency range.



The ERKA.phon bell is available in the ADULT and CHILD sizes. A ring as protection against cold is available for the membrane side as an option.



The bell – available in brass/ chrome-plated brass with high-gloss finish/chrome-plated with matt finish or anodized aluminum – is quick and easy to clean.



The tube is available in many color variations.



To ensure that signals captured by the bell are transmitted to the ears without interference, the spring was not integrated in the tube duct.



The ERKA.phon is also available in a special version for teaching applications.



	Material · Surface	black	grey	blue	green	burgundy	red	gright green	violet	dark grey
Type: double Size: adult	Brass · bright chrome	541.00010	541.00000	541.00020	541.00055	541.00060	541.00030	541.00051	541.00080	541.00005
	Brass · matt chrome	541.01010	541.01000	541.01020	541.01055	541.01060	541.01030	541.01051	541.01080	541.01005
	Aluminium · anodized	542.00010	542.00000	542.00020	542.00055	542.00060	542.00030	542.00051	542.00080	542.00005
Type: double Size: child	Brass · bright chrome	545.00010	545.00000	545.00020	545.00055	545.00060	545.00030	545.00051	545.00080	545.00005
	Brass · matt chrome	545.01010	545.01000	545.01020	545.01055	545.01060	545.01030	545.01051	545.01080	545.01005
Type: flat Size: adult	Brass · bright chrome	543.00010	543.00000	543.00020	543.00055	543.00060	543.00030	543.00051	543.00080	543.00005
	Brass · matt chrome	543.01010	543.01000	543.01020	543.01055	543.01060	543.01030	543.01051	543.01080	543.01005
	Aluminium · anodized	544.00010	544.00000	544.00020	544.00055	544.00060	544.00030	544.00051	544.00080	544.00005
Type: flat Size: child	Brass · bright chrome	547.00010	547.00000	547.00020	547.00055	547.00060	547.00030	547.00051	547.00080	547.00005
	Brass · matt chrome	547.01010	547.01000	547.01020	547.01055	547.01060	547.01030	547.01051	547.01080	547.01005

The **ERKA.phon** is available in different versions:

one for using in a mr-room and one special version for teaching applications; the matching tubes are available in many different colors, for example in special coloration pink  or light blue 

Spare-part overview

[illegible]

ERKA. Pressure Infusion Cuffs. **The latest in transfusion and infusion.**

We designed our **ERKA.** pressure infusion cuffs for daily use in the operating theatre, intensive care and rescue vehicles.

The use of components in the **ERKA.** SWITCH 2.0 product family with all its benefits in terms of handling and precision is a new feature. **ERKA.** also offers groundbreaking innovations with the new line of disposables – typical for **ERKA.**

High Pressure Cuffs (HPC)



ERKA. SWITCH 2.0 HPC

Efficient pressure application.

The reusable pressure infusion cuff from the innovative **ERKA.** Switch 2.0 product family features sophisticated details and efficient handling.

The outer cover is washable, the mesh fabric is extremely durable and resilient, the bulb is made of high-quality PU and all rubber components are naturally latex-free.



Available for use with 500 ml and 1000 ml transfusion/infusion bags.




The requirements for use in an MR room are met.



The polycarbonate material is extremely easy to clean as well as impact-resistant based on internationally recognised standards such as ANSI/AAMI SP10. The scale with a diameter of 56 mm is optimised for legibility.



Ergonomic and balanced, made of lightweight materials and equipped with a sensitive threaded discharge valve, the **ERKA.** SWITCH 2.0 HPC guarantees maximum ease of use.



Replacing the load cell? Child's play. Exchange it for a new, calibrated load cell in a few simple steps. Without the help of qualified medical-technical personnel!

	Cuff · Type · Color	Size 1: 10-15 cm	Size 2: 14-21 cm	Size 3: 20,5-28 cm	Size 4: 27-35 cm	Size 5: 34-43 cm	Size 6: 42-54 cm
SWITCH 2.0 Simplex, one cuff	superb · D-Ring · grey	-	-	294.46891	294.40891	294.44891	-
SWITCH 2.0 Comfort, one cuff	superb · Rapid · grey	293.28493	293.27493	293.26493	293.20493	293.24493	293.29493
	superb · Hook · grey	-	-	293.26293	293.20293	293.24293	-
	superb · D-Ring · grey	-	-	293.46893	293.40893	293.44893	-
	smart · Rapid · grey	293.38393	293.37393	293.36393	293.30393	293.34393	293.39393
SWITCH 2.0 Comfort with cuff-set 1: Size 1-3	superb · Rapid · grey	293.71393	Part of the set	Part of the set	-	-	-
	smart · Rapid · grey	293.74393	Part of the set	Part of the set	-	-	-
SWITCH 2.0 Comfort with cuff-set 2: Size 3-5	superb · Rapid · grey	-	-	293.71293	Part of the set	Part of the set	-
	superb · Hook · grey	-	-	293.72293	Part of the set	Part of the set	-
	superb · D-Ring · grey	-	-	293.74293	Part of the set	Part of the set	-
	smart · Rapid · grey	-	-	293.73293	Part of the set	Part of the set	-
SWITCH 2.0 Comfort with cuff-set 3: Size 1-6	superb · Rapid · grey	293.71193	Part of the set	Part of the set	Part of the set	Part of the set	Part of the set
	smart · Rapid · grey	293.73193	Part of the set	Part of the set	Part of the set	Part of the set	Part of the set
	Size 1: 500 ml	Size 2: 1000 ml					
SWITCH 2.0 HPC Pressure infusion cuff	822.00000	832.00000					

Accessories for SWITCH 2.0

	Housing for Comfort	Housing for HPC and Simplex (with hook)	Hook (separat)	safety clip (separat)	Body with ball for Comfort
SWITCH 2.0	000.30200	000.30201	995.00139	995.00140	000.30100

ERKA. blood pressure cuffs.
Universal, comfortable and green.

The virtually limitless selection of **ERKA.** blood pressure cuffs – distinguished by unrivalled patient comfort and hygienic, ecologically convincing materials – is not only used in our own products. Our OEM partners worldwide in the fields of NIBP monitoring, 24-hour blood pressure measurement and ergometry have been putting their trust in us as a supplier of custom blood pressure measurement cuffs for decades.



Cuffs

ERKA. Green Cuff - Reusables

Superb, smart and soft.

Superb because all ERKA. Green Cuff products are PVC and latex-free, every one is individually tested prior to delivery and we guarantee 50,000 inflation cycles.

Smart because all ERKA. Green Cuff products are equipped with clearly marked size and positioning features for the most accurate possible measurements.

Soft because the supple, easy to disinfect, biocompatible materials ensure maximum patient comfort.



The ERKA. Green Cuff is available in 3 versions: SUPERB (reusable, consisting of cover and bulb), SMART (reusable, cover with integrated bulb) and SOFT (disposable, cover with integrated bulb).



The ERKA. Green Cuff Superb is available in the RAPID, HOOK and BRACKET functional versions. Smart and Soft cuffs only in the RAPID version.



ERKA. Green Cuff series products are available in two standard colours.



Thanks to the implementation of the international ANSI/AAMI SP10 standard, ERKA offers blood pressure measurement cuffs in 6 sizes: 10-15 cm | 14-21.5 cm | 20.5-28 cm | 27-35 cm | 34-43 cm | 42-54 cm



DISPOSABLES are available in ERKA. Green Cuff Soft with a fleece side and a side that can be wiped, and in ERKA Green Cuff Soft² with two fleece sides.

	Cuff · Type · Color	Size 1: 10-15 cm	Size 2: 14-21 cm	Size 3: 20,5-28 cm	Size 4: 27-35 cm	Size 5: 34-43 cm	Size 6: 42-54 cm
Single-channel-tube cuff for - SWITCH 2.0 Comfort - Kobold	superb · Rapid · green	025.28483	025.27483	025.26483	025.20483	025.24483	025.29483
	superb · Rapid · grey	025.28493	025.27493	025.26493	025.20493	025.24493	025.29493
	superb · Hook · green	-	-	025.26283	025.20283	025.24283	-
	superb · Hook · grey	-	-	025.26293	025.20293	025.24293	-
	superb · D-Ring · green	-	-	025.46883	025.40883	025.44883	-
	superb · D-Ring · grey	-	-	025.46893	025.40893	025.44893	-
	smart · Rapid · grey	025.38393	025.37393	025.36393	025.30393	025.34393	025.39393
Double-channel-tube cuff for - Perfect Aneroid - Profi - SWITCH 2.0 Simplex - Vario - Klinik ERKA .meter - ERKA .meter 3000	smart · Rapid · grey	025.38391	025.37391	025.36391	025.30391	025.34391	025.39391
	superb · Hook · grün	-	-	025.26281	025.20281	025.24281	-
	superb · Hook · grey	-	-	025.26291	025.20291	025.24291	-
	superb · D-Ring · grün	-	-	025.46881	025.40881	025.44881	-
	superb · D-Ring · grey	-	-	025.46891	025.40891	025.44891	-
	superb · Rapid · green	025.28481	025.27481	025.26481	025.20481	025.24481	025.29481
	superb · Rapid · grey	025.28491	025.27491	025.26491	025.20491	025.24491	025.29491

Please find our cuff accessories and tubes at the spare parts.

The **ERKA** Green Cuff:

Durable: The bladders for our new cuffs consist of a patented German PU material that is extremely durable and resilient.

Hygienic: The cover of the new **ERKA**. cuffs is excellent to clean and to disinfect. This allows a maximum of hygiene. All of our new cuffs are free from latex, PVC and DEHP and therefore skin compatible and of course biocompatible.

Practicable: The new design is user friendly and easy to handle in daily use. Our assortment comprises rapid, hook and bracket cuffs in up to 6 different sizes and two colors (grey and green).



That easy you change the new **ERKA** Green Cuff:

1st: Pull off the tube at the bladders exit of the old cuff. This one stays with the device and is due to its persistence reusable.

2nd: Fix the tube to the new replacement-cuff. This one always features the **ERKA**. standard connector. Then heat the end of the tube in a water bath to facilitate the matching.

These instructions are valid for both one-tube and two-tube versions.

ERKA. Green Cuff – Reusables | Cuff-Sets

Cuff-Sets can only be used with **ERKA.** devices equipped with the **ERKA.** quick-fastener.

	Cuff · Type · Color	Size 1: 10-15 cm	Size 2: 14-21 cm	Size 3: 20,5-28 cm	Size 4: 27-35 cm	Size 5: 34-43 cm	Size 6: 42-54 cm
Single-channel-tube cuff for - SWITCH 2.0 Comfort - Kobold	superb · Rapid · grey	024.28493	024.27493	024.26493	024.20493	024.24493	024.29493
	superb · Hook · grey	-	-	024.26293	024.20293	024.24293	-
	superb · D-Ring · grey	-	-	024.46893	024.40893	024.44893	-
	smart · Rapid · grey	024.38393	024.37393	024.36393	024.30393	024.34393	024.39393
Double-channel-tube cuff for - Vario - Klinik ERKA. meter	superb · Hook · green	-	-	024.26281	024.20281	024.24281	-
	superb · D-Ring · green	-	-	024.46881	024.40881	024.44881	-
	superb · Rapid · green	024.28481	024.27481	024.26481	024.20481	024.24481	024.29481
	smart · Rapid · grey	024.38391	024.37391	024.36391	024.30391	024.34391	024.39391

Please find our cuff accessories and tubes at the spare parts.

Set cuffs:

The most frequent mistake made in blood pressure measuring is the choice of a cuff that is too small or too large. The right size of the cuff is guaranteed with our **ERKA.** set cuff:

Being equipped with a colored tube – that is conform to international standards – and the **ERKA.** standard connector, you can change the cuffs in a matter of seconds and also recognize the right size of the cuff at a glance.

Of course you can exchange each **ERKA.** product with a quick fastener. Products, which you have bought with a cuff set, feature the **ERKA.** quick fastener already.



The **ERKA.** quick fastener
(see picture) as spare part:
995.0001

The **ERKA.** quick fastener permits the exchange
of the cuffs in a matter of seconds.

ERKA. Classic Cuff - Reusables

The classic in green.

The reusable ERKA. Classic Cuff features a rubber bulb and cotton or PU-coated cotton cover. It is available in the rapid, hook or bracket functional versions, covering all common requirements with sizes from newborn to thigh.

Thanks to their characteristics, the materials used guarantee easy cleaning, ultimate surface hygiene and maximum wearer comfort.



	Cuff	Type	Color	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6
Single-channel tube version - Kobold	Cotton	· Rapid	· green	025.08303	025.07303	025.06303	025.00303	025.94303	025.04303
	wipeable	· Rapid	· green	025.08403	025.07403	025.06403	025.00403	025.94403	025.04403
	Cotton	· Hook	· green	-	-	025.06103	025.00103	025.94103	025.04103
	wipeable	· Hook	· green	-	-	025.06203	025.00203	025.94203	025.04203
Double-channel tube version - Vario - Klinik ERKA .meter - ERKA .meter 3000	Cotton	· Rapid	· green	025.08301	025.07301	025.06301	025.00301	025.94301	025.04301
	Cotton	· Hook	· green	025.08401	025.07401	025.06401	025.00401	025.94401	025.04401
	wipeable	· Rapid	· green	-	-	025.06101	025.00101	025.94101	025.04101
	wipeable	· Hook	· green	-	-	025.06201	025.00201	025.94201	025.04201
Double-channel tube version - Perfect Aneroid - Profi	Cotton	· Rapid	· green	025.08301	025.07301	025.06302	025.00302	025.94302	025.04302
	wipeable	· Rapid	· green	025.08401	025.07401	025.06402	025.00402	025.94402	025.04402
	Cotton	· Hook	· green	-	-	025.06102	025.00102	025.94102	025.04102
	wipeable	· Hook	· green	-	-	025.06202	025.00202	025.94202	025.04202

**ERKA. Made in Germany. Material of the highest quality.
Ultimate quality. Absolute perfection.**

As the blood pressure measurement pioneer, our products have to meet the highest standards. The name **ERKA** has represented the pinnacle of reliability and professionalism for 120 years.

For example with the "Perfect Aneroid". It has to pass the **ERKA** Quality Management Check - QMC for short - before it is delivered. A total of 45 individual inspections, from the material inspection the comprehensive verification of functionality, ensure perfection you can rely on.

Nevertheless, if you are looking for a replacement part for one of our products, we are always there for you.



Spare Parts

	Material · Color	Perfect Aneroid	Profi	Kobold	Switch 2.0 Comfort / Simplex	Vario	E3000/eco	Klinik ERKA.meter
Case	green	000.14000	000.14000	000.14000	-	-	-	-
	black	000.14200	000.14200	000.14200	-	-	-	-
	white	-	-	-	000.14404	-	-	-
Quick fastener connector	Brass chrome-plated	-	-	995.00001	995.00001	995.00001	-	995.00001
Ball	PVC · green	000.20200	-	-	-	000.20200	000.20200	000.20200
	PVC · black	-	000.20400	000.20400	-	000.20400	000.20400	000.20400
	PVC · grey	-	-	-	000.20900	000.20900	000.20900	000.20900
Milled screw	Brass chrome-plated	000.35200	000.35200	000.35200	-	000.35200	000.35200	000.35200
Hygienic protection role		000.04100	000.04100	000.04100	000.04100	000.04100	000.04100	000.04100
Pressure reducing valve	Brass chrome-plated	-	-	-	-	000.51000	000.51000	000.51000
	PC-Polyman	-	-	-	-	000.21000	000.21000	000.21000
Valve + ball green	Brass	-	-	-	-	000.20500	000.20500	000.20500
Valve + ball black	Brass	-	-	-	-	000.20700	000.20700	000.20700
Valve + ball grey	PC-Polyman	-	-	-	-	000.21100	000.21100	000.21100
Extension tube (per meter, single tube)	Neopren · green	000.58400	000.58400	000.58400	000.58400	000.58400	000.58400	000.58400
	Silicon · grey	-	-	995.00053	995.00053	995.00053	-	995.00053

	Material · Farbe	Perfect Aneroid	Profi	Kobold	Switch 2.0 Comfort / Simplex	Vario	E3000/eco	Klinik ERKA .meter
Connector male	Brass chrome-plated	-	-	-	-	000.75502	000.75502	000.75502
Connector female	Brass chrome-plated	-	-	-	-	000.75501	000.75501	000.75501
Connector complete	Brass chrome-plated	-	-	-	-	000.75500	000.75500	000.75500
End cap (ME = 10 Pieces)	Brass chrome-plated	000.00500	000.00500	000.00500	-	000.00500	000.00500	000.00500

Cuff equipement and tubes

	Perfect Aneroid	Profi	Kobold	Switch 2.0 Comfort	Switch 2.0 Simplex	Vario	E-3000/eco	KlinikERKA meter
single tube · grey	-	-	995.00053	995.00053	-	-	-	-
double tube · green	-	-	-	-	-	995.00070	995.00070	995.00070
double tube · grey	-	-	-	-	on request	on request	on request	on request
twin tube · green	000.58407	000.58407	-	-	-	-	-	-
twin tube · grey	000.58408	000.58408	-	-	-	-	-	-
spiral coiled tube (Silicon) · green	-	-	-	-	-	000.56000	-	000.56000
spiral coiled tube (Silicon) · grey	-	-	-	-	-	000.56200	-	000.56200

[illegible]

General Conditions of Sale and Delivery

1. Acceptance of the Conditions of Sale and Delivery

„Our conditions of delivery and payment are exclusively applicable and agreed to by our customer when placing an order and also apply to future transactions, even if they are not exclusively referred to, but received by the customer when we confirmed the receipt of an order. If the order is placed deviating from our conditions, our delivery and payment conditions still apply, even if we do not raise objections. So deviations are only applicable if they are expressly recognized by us in writing.“
„We are entitled to assign the claims based on our business connection.“

2. Placing of Orders

The offers of the supplier are not binding. The contract is concluded when the supplier has confirmed the acceptance of the order in writing within a period of 4 weeks or when the delivery has been effected. All agreements have to be made in writing; this also applies to supplements, amendments and collateral agreements. The customer is liable for the correctness of the documents to be supplied by him, especially samples and drawings. The details, drawings, illustrations and specifications of services of the supplier included in the catalogues, price lists or documents being part of the offer are approximate values customary in the trade. Quality, design and colour are subject to customary and material deviations, if this is unavoidable due to the situation on the raw material market or due to technical reasons.

3. Delivery Delivery

dates or periods, which can be stipulated with or without binding force, have to be indicated in writing. Delivery periods start on conclusion of the contract. Delivery deadlines are deemed to be met if the goods have left the factory/warehouse at the agreed date or at the date when the goods are ready for dispatch, as far as the dispatch or the delivery of the goods is delayed by circumstances which are under the control of the customer. In case of subsequent amendments of the contract by the customer which influence the delivery time, it is extended appropriately. In case of delay in delivery, an appropriate deadline for delayed delivery has to be set – starting on the date of receipt of the written notice of default made out by the customer or, in the case of delivery periods set by calendar, when these expire. This also applies to partial deliveries. Disturbances of our business activities, especially strikes, lockouts, lack of manpower, even on the part of our subcontractors and manufacturing firms, lacking dispatch possibilities and procurement of raw material, cases of force majeure or cases, in which we are not supplied by our subcontractors and the non-supply has not negligently been caused by us, change the delivery periods by the duration of the interference with the performance of the obligation caused by these circumstances or discharge us from the delivery commitment to the exclusion of any claims for damages as far as the removal of the obstacle is not foreseeable with respect to time. The existence of such circumstances and obstacles which are beyond our control has to be proved by us. If the abovementioned obstacles occur on the part of the customer, the same legal consequences also apply to his obligation to take delivery. Supplier as well as customer are obliged to inform the other party of the beginning and the end of the above-mentioned kind of obstacles without delay.

4. Pricing

The prices are understood ex works plus value-added tax and do not include packing and postage. If prices are not listed or only with the reservation „current list price“, the list prices valid on the delivery day will be invoiced. This only applies to delivery periods in excess of 4 months and to price adjustments up to 4%. Higher price adjustments require new price agreements. If such an agreement does not exist, the customer is entitled to withdraw from the contract. In case of free delivery it is up to the supplier to determine the mode of dispatch. If goods are sent by fast freight or express, the customer bears the difference between the charges for freight and those for the accelerated dispatch, even in case of free delivery.

5. Conditions of Payment

All invoices are payable within 30 days after invoicing without discount. For payments within 8 days, a 2% discount may be deducted, provided that there is no delay of payment relating to previous deliveries. When the credit period is exceeded, the supplier is entitled to charge default interest amounting to 2% above the corresponding discount rate of the Deutsche Bundesbank. It is to be estimated at a higher or lower rate, if the supplier furnishes proof of a higher interest load or the customer of a lower load. Money orders, cheques and bills of exchange are only accepted after a special agreement and only on account of payment, invoicing all collection and discounting charges. If, after the conclusion of the contract, the supplier learns about facts considerably deteriorating the financial conditions of the customer and suited to endanger his claim to counterperformance according to his best and commercial judgement, he can demand furnishing of a suitable security within an appropriate period of time or performance in return for counter-performance up to the performance date. If the customer fails to comply with the justified demand of the supplier or does not comply with it in time, the supplier is entitled to withdraw from the contract or to claim damages for non-performance. If the customer is in arrears with a part performance, the supplier can accelerate the whole residual claim immediately. In case of delay of performance resulting from a considerable deterioration of the financial situation, the supplier can withdraw from the contract or claim damages for non-performance without granting a grace period. If the delay of performance is not property-related, the supplier can demand withdrawal from the contract or damages for non-performance after fruitless expiry of an adequately set period of time. The customer can only set off against claims of the supplier or effectively enforce a lien, if his counterclaim has become res judicata or has been uncontested or if the right of retention results from the contractual relationship which these conditions are the basis of.

6. Dispatch and Passing of Risk

The goods are dispatched ex factory, if not otherwise provided in the agreement, without the cheapest mode of dispatch being obligatory. Even in case of freight-free delivery, the risk is passed to the customer when the goods are handed over to the person entrusted with the dispatch or when being loaded onto a vehicle of the supplier. The supplier is entitled, but not obliged, to insure deliveries on behalf and for account of the customer. If goods are damaged in transfer, the customer has to cause the competent authorities to ascertain the facts without delay and to inform the supplier. Partial deliveries on a reasonable scale are permissible.

7. Reservation of Title

The supplier reserves his title to the delivered goods until all claims out of the business relation are satisfied, including the claims accruing in future even from contracts concluded with the customer at the same time or later on. The customer is entitled to sell the goods in the ordinary course of business as long as he fulfils the obligations arising out of the business relations with the supplier in good time. However, he is neither allowed to pledge the reserved goods nor assign them by way of security. He is obliged to secure the rights of the supplier when reselling the reserved goods on credit. If the customer gets into arrears or if duties are neglected resulting from the reservation of title, the supplier is entitled to demand the temporary return of the reserved goods at the expense of the customer and to collect them from the customer even without exercising the right to rescind the contract and without granting a grace period. The customer has not the right of possession in that case. Possible processing or manufacturing of the reserved goods is always done by the customer for the supplier. If the reserved goods are processed or inseparably mixed with other objects not belonging to the supplier, the supplier acquires ownership of the new object in proportion of the invoice value of the reserved goods to the other processed or mixed objects at the time of the processing or mixing. If goods of the supplier are merged or inseparably mixed with other movable objects to form a homogeneous object and if the other object is considered to be the main object, it is regarded as agreed that the customer proportionally transfers co-ownership to the supplier, as far as the main object belongs to him. All claims and rights arising from the sale of goods, of which the supplier is entitled to an ownership, are assigned by the customer to the supplier for security already. The supplier herewith accepts the assignment. If the customer is in arrears or if other obligations are infringed arising from

the agreed reservation of title, the supplier is entitled to disclose to the purchasers of the customer the assignment of the claim of the customer to the supplier and to collect the claims. The customer has to inform the supplier without delay about execution measures of third parties levied upon the reserved goods or upon the claims assigned to the supplier or other securities, handing over the documents necessary for an intervention. This also applies to prejudices of another kind. The supplier obliges himself to release the securities he is entitled to according to the above provisions upon request of the customer in so far as the value of the goods assigned for security exceeds the claims to be secured by more than 20%.

8. Warranty

If the delivered goods are defective or if warranted characteristics are missing, the supplier must exclusively warrant the goods within 24 months after the day of delivery in such a way that he either remedies the defect free of charge or subsequently delivers goods free of defects. Other claims of the customer due to defects or the absence of warranted characteristics even for damages, especially for consequential harms caused by a defect, are excluded. If the customer is a trader, he has to give notice of obvious defects immediately. If he is a non-trader, the period of giving notice of a defect is 10 days after taking delivery of the goods. If the customer fails to observe these time limits, claims for damages cannot be asserted any more. Notice of non-obvious defects must be given immediately after the defect becomes noticeable. Notices of defects have to be made in writing. For observing the time limit it is sufficient to post the notice of defect in time. The requirements concerning form and time limit for notices of defect also apply to effected partial deliveries. For remedying all defects or for substitute deliveries which appear necessary according to the reasonably exercised discretion of the supplier, the customer has to grant to the supplier after informing him, the necessary time and opportunity, otherwise the supplier will be relieved of the liability for defects. Only in urgent cases of endangering the operational safety or for warding-off disproportionately severe damages, in which the supplier has to be informed without delay, or if the supplier is in default of remedying the defect, the customer is entitled to remedy the defect himself or to have it remedied by third persons and to demand reimbursement of the necessary costs from the supplier. If the supplier lets expire the reasonable grace period granted to him for remedying a defect or for making a substitute delivery, or if remedying the defect or making a substitute delivery has failed twice, the customer is entitled to withdraw from the contract or to demand an abatement of the purchase price. The assertion of claims for damages, especially due to late substitute delivery or failure to make a substitute delivery, late or failed remedy of defects, especially for consequential damages, is excluded. For remedying of defects and substitute deliveries the supplier is liable to the same extent as for the original delivery item. For substitute deliveries the original warranty period applies. The supplier is not obliged to remedy a defect or to make a substitute delivery, as long as the customer does not fulfil his obligations under the contract.

9. Place of Performance and Jurisdiction

The head office of the supplier is the place of performance for all obligations arising out of the contractual relationship. For all legal actions, even for special procedures deciding claims arising out of a bill of exchange and for actions for assertion of a claim concerning payment of a cheque, the exclusive „place of jurisdiction“ is the „domicile of the company“ if the customer is a registered trader, a legal person under public law or of the special revenue fund under public law. This place of jurisdiction is equally applicable if the customer has no general place of jurisdiction within the country, transfers his domicile or habitual residence after conclusion of the contract outside the country or if his domicile or habitual residence is not known when legal proceedings are instituted. The contract is subject to German law as far as it has been concluded with a foreign party. The Uniform Law on the International Sale of Goods as well as the Uniform Law on the Formation of Contracts for the International Sale of Goods are not applicable. These conditions supersede previous conditions of sale and delivery. If one of the provisions of these conditions of sale and delivery should become invalid, this shall have no effect upon the other conditions. Personal data arising out of our terms and conditions of business will be stored in our EDP.

ERKA.
THE ORIGINAL
Made in Germany since 1889