


VIGMED® CLiP® PORTED

SAFE & *consistent*

Developed with users and patients in mind, CLiP products are configured to give excellent penetration characteristics. Patients and their conditions may vary, but our products must remain consistent to give users the confidence that safe first stick success is achievable, every time.



PRECISE TIP CONFIGURATION

The consistent interplay between the needle and catheter tip is based on years of development experience to give ideal penetration characteristics.

ULTRA SHARP, ULTRA POLISHED

To help users achieve consistent penetration, all needles are of superior Japanese design. They are electropolished and back-cut, making them extremely sharp. They cut rather than force a hole, which helps healing and leaves minimal scarring.

AUTOMATIC PASSIVE SAFETY

The safety mechanism is automatic so no additional action is required by the user during insertion. There is a precisely designed bump on the needle, which ensures that the safety mechanism is activated.

AVAILABLE IN PUR AND FEP

Our catheters are made from medical grade virgin polymers. You have the option between PUR and FEP. PUR becomes softer in contact with moisture and heat, which in some situations is considered better for patient care. FEP is like Teflon and therefore virtually friction-free.

PRODUCT CODE	DESCRIPTION	MODEL	GAUGE	ØxLENGTH mm	COLOR CODE	FLOW RATE ml/min	PACK SIZE pcs
VP14 45 11	CLiP Ported Safety I.V. Catheter	PUR	14	2,0 x 45	Orange	290	50
VP16 45 11	CLiP Ported Safety I.V. Catheter	PUR	16	1,7 x 45	Grey	200	50
VP17 45 11	CLiP Ported Safety I.V. Catheter	PUR	17	1,5 x 45	White	140	50
VP18 45 11	CLiP Ported Safety I.V. Catheter	PUR	18	1,2 x 45	Green	100	50
VP18 32 11	CLiP Ported Safety I.V. Catheter	PUR	18	1,2 x 32	Green	110	50
VP20 32 11	CLiP Ported Safety I.V. Catheter	PUR	20	1,0 x 32	Pink	64	50
VP22 25 11	CLiP Ported Safety I.V. Catheter	PUR	22	0,8 x 25	Blue	38	50
VP24 19 11	CLiP Ported Safety I.V. Catheter*	PUR	24	0,7 x 19	Yellow	22	50
VP14 45 01	CLiP Ported Safety I.V. Catheter	FEP	14	2,0 x 45	Orange	290	50
VP16 45 01	CLiP Ported Safety I.V. Catheter	FEP	16	1,7 x 45	Grey	200	50
VP17 45 01	CLiP Ported Safety I.V. Catheter	FEP	17	1,5 x 45	White	140	50
VP18 45 01	CLiP Ported Safety I.V. Catheter	FEP	18	1,2 x 45	Green	100	50
VP18 32 01	CLiP Ported Safety I.V. Catheter	FEP	18	1,2 x 32	Green	110	50
VP20 32 01	CLiP Ported Safety I.V. Catheter	FEP	20	1,0 x 32	Pink	64	50
VP22 25 01	CLiP Ported Safety I.V. Catheter	FEP	22	0,8 x 25	Blue	38	50
VP24 19 01	CLiP Ported Safety I.V. Catheter*	FEP	24	0,7 x 19	Yellow	22	50

*With needle notch

MATERIALS

Needle: Stainless steel | Catheter: PUR, FEP | Other parts: PP, PE, LDPE, HDPE, POM | Silicone | Silicone oil