



Woburn Ultimate Profiling Bed

Overview

Product Name: Woburn Ultimate Profiling Bed

Product Code: HLB798/12

Description: The Woburn Ultimate Profiling Bed is our higher capacity bed within the Woburn range. Featuring triple aluminium side rails and under-bed strengthening, the Ultimate excels in patient support and ease of infection control. Available in 1200mm wide, with a 60 Stone weight capacity, this bed is able to support most bariatric patients. Comfort adjustments are made with ease due to the powerful lifting components and the high quality electronic control system.

Safe Working Load: 380kg / 60 Stone

Max. Patient Weight: 300kg

Overall Dimensions: 2195 x 1340mm (1200mm mattress base)

Max. Mattress Depth: 260mm

Adjustable Height: 275 - 690mm

Angle of Back Rest: Up to 70 degrees

Angle of Knee Break: Up to 30 degrees

Bed Weight: 157kg + 16kg with rails

Key Information:

- 4 Section Profiling
- Electronically Powered Back Rest & Knee Break
- Manual Leg Rest Adjustment
- Wraparound Wooden Head and Foot Boards
- Aluminium Side Rails
- Auto-Regression
- Trendelenburg Function - Optional
- 200mm Extension - Optional
- Lockable Handset
- Low Voltage – 29v; 9v Battery Backup
- Comprehensive Manual
- CE marked
- Warranty: 5yrs

Medguard Professional Health Care

Unit 6/7, Block 4, Ashbourne Business Park, Ashbourne, Co Meath

Tel: + 353 (0) 1 835 2411

Email: service@medguard.ie

Fax: +353 (0) 1 969 5009

Website: www.medguard.ie