

**Product Code: CSNT/4**

**Description: Case Notes Trolley, Eight Removable Shelves**



### Features

- Manufactured from mild steel
- Secured with CAM lock
- Powder coated blue
- Eight removable shelves
- Lower storage section
- Mounted on 4x 100mm swivel Tente ® Castors
- Protective rubber buffers
- White MDF top

### Dimensions (WxDxH)

610mm X 461mm

HPRA Wholesalers authorisation No. W00426/00001

Tel: +353 (0) 1 835 2411

Email: [sales@medguard.ie](mailto:sales@medguard.ie)

Unit 28B Ashbourne Business Centre, Ashbourne, Co Meath, A84 WA49