

Product Data Sheet: 112055 / 550mm

## Nymapro Trombone lift and lock steel hinged grab rail



### Product Detail:


This hinged support rail has been designed with a lift and lock folding mechanism and a integral toilet roll holder.

#### Features include:

- Lift and lock folding mechanism
- Integral toilet roll holder
- Robust steel construction
- Durable epoxy polyester coating
- Doc M compliant (*800mm ONLY*)

MAX RATED LOAD: 120 kg

Reference (specify): 112055 / 550mm



Order Code	Size (A)
112055	550mm
DDGR-B	800mm

PLEASE NOTE: Every effort has been made to ensure as much accuracy as possible. However, due to print restrictions exact colour match is not always possible. Featured products may also be changed from time to time, but are correct at time of publication.

\* Terms and Conditions apply - [www.nymas.co.uk/warranties](http://www.nymas.co.uk/warranties)