



Hausted Video Imaging Chair

A Perfect Solution when Performing
Modified Barium Swallow Examination

The Video Imaging Chair is Specially Designed
for Safe Patient Transport whilst undergoing
Fluoroscopic Examination

Operational Benefits

- Specially Design of Back Section accommodates the limited space available when using Fluoroscopic Equipment
- Radiolucent Back Section provides a clear unobstructed Radiological Image
- Back Section Transport Attachment - Quick and Easy Removal
- Full Range Patient Positioning is achieved using the Backrest and Leg Sections, adjust from 90° Upright to Full Supine, and 10° Trendelenburg
- The Rotating Seat Section enables quick repositioning from AP to Lateral Imaging View
- Stainless Steel Base Cover - Seamless and Easy Clean
- Four Wheel Brake & Steer System provides Triple Mode Operation: Steer / Neutral / Brake



Optional Accessories

- Chair Foot Extension
- Patient Tray
- Paediatric Patient Positioners - Large
- Paediatric Patient Positioners - Medium
- Paediatric Patient Positioners - Small
- Patient Tray
- Vertical Oxygen Tank Holder
- Arm Board
- Paper Roll Holder
- Mayo Stand with Tray



(Above) AP Projection

(Left) Chair Rotated by 90° Lateral Projection

Patient Outcome & Safety Benefits

- VIC can raise in the supine position to safely transfer a patient from Stretcher or Bed to Chair
- Reduced risk of injury to transporting personnel
- High-density, non-conductive vinyl-covered foam chair pads support the patient comfortably. Padded Arm Rests for comfort
- Pneumatic Assisted Backrest Lever adjusts both sides of the chair from supine to 90 degrees for patient comfort and procedure positioning. An additional lever controls a 10 degree Trendelenburg
- Leg Section and Fold-Down Footrest adjust, maximising Lower Body Comfort and Stability

Specifications

- Seat Cushion Width 56 cm
- Overall Length (Supine) 195 cm
- Overall Length (Upright) 106 cm
- Height Adjustment 60 - 80 cm
- Wheel Base 63 x 57 cm
- Pneumatic-Assist Back-rest 0 - 90°
- Trendelenburg 10°
- Maximum Weight Capacity 147.5 kg

Medguard Professional Healthcare

Tel: +353 (0) 1 835 2411

Email: sales@medguard.ie

Fax: 01 969 5009

Website: www.medguard.ie

Unit 6/7, Block 4, Ashbourne Business Park, Ashbourne, Co. Meath, A84 VK52