



## Specification Sheet - Advancis Medical Wound Probes (50)

### CR4236 - Advancis Medical Wound Probes (50)



Description	Advancis Medical Wound Probe
Length	10cm
Sterile	Yes
Use	Wound assessment stick designed to determine the depth of a wound
Indications	Pressure ulcers, sinuses, wounds with an unknown depth
Contra-indications	Yes
Use	Wound assessment stick designed to determine the depth of a wound

#### Directions for use:

Use the probe to assess the wound by inserting into the wound. If the wound opening is smaller than the width of the probe a different form of assessment should be used and the probe should not be inserted.

Wound Probe® is a single use product. If used on more than one patient, cross-contamination or infection may result.

Opening the dressing pack compromises the sterile barrier. Any unused dressings should be discarded.