



BLUE TRANSFER BOARD

The safe way to transfer patients

A legal requirement in some countries, taken for granted in others – working conditions should be as comfortable as possible. Healthcare workers should avoid lifting when transferring patients. **SAMARIT Products** assist with transfer, without the need to lift. The products are designed to facilitate the transfer without using a lot of effort, thus making the work easier for the carers and kinder to their backs. Can be ordered as OEM product.

BLUE LONG WIDE
BLUE SHORT NARROW
BLUE LONG NARROW
BLUE SHORT WIDE



Medguard Professional Healthcare

Tel: +353 (0) 1 835 2411

Email: sales@medguard.ie

Fax: 01 969 5009

Website: www.medguard.ie

Unit 6/7, Block 4, Ashbourne Business Park, Ashbourne, Co. Meath, A84 VK52

SAMARIT

Samarit transfer systems

BLUE LONG WIDE

foldable



Suitable for radiology, operating theatre, intensive care units and all special wards.

Size: 180 x 50 cm / 89 x 50 cm, 2.9 kg
70 x 20 inch / 35 x 20 inch 6.4 lbs

SWL*: 240 kg / 529 lbs

Order No. 440.8300
UDI/GTIN: 7640144460228

BLUE SHORT NARROW

non-foldable



For ambulance, operating theatres, wake-up stations and emergencies.

Size: 100 x 42 cm, 1.4 kg
39 x 17 inch 3.1 lbs

SWL*: 150 kg / 331 lbs

Order No. 440.8301
UDI/GTIN: 7640144460235

BLUE LONG NARROW

foldable



Suitable for operating theatre, X-Ray, various wards.

Size: 180 x 42 cm / 89 x 42 cm, 2.1 kg
70 x 17 inch / 35 x 17 inch 4.6 lbs

SWL*: 180 kg / 397 lbs

Order No. 440.8305
UDI/GTIN: 7640144460259

*SWL = Safe Weight Limit

4 different models to suit every hospital application

- High performance
- Radioluscent
- Non-toxic
- Anti-bacterially surface
- Safe and comfortable for patients
- Prevents back injuries to care workers
- Skin-friendly (no latex)
- Easy to clean and disinfect – hygiene approved
- Disposable covers available
- Tear resistant nylon cloth
- Bridges gaps up to 20 cm
- Lightweight, easy suspendable

BLUE SHORT WIDE

non-foldable



For urology operating theatres, gynaecology operating theatres, wake-up stations and emergencies.

Size: 80 x 50 cm, 1.4 kg
3.15 x 2 inch 3.1 lbs

SWL*: 150 kg / 331 lbs

Order No. 440.8304
UDI/GTIN: 7640144460242