

WOBURN ULTRA-LOW BED



Product Code:

HLB799

Description:

- A fully functioning four- section Ultra-Low profiling bed.
- Platform can be lowered to 67mm off the floor for patient safety.
- Electronically setting of height, back rest or knee brake is simple, using the lockable handset.
- The Woburn Ultra-Low has a Trendelenburg function as standard.
- Optional bed levels, side rails, lighting and accessories available

Safe Working Load:

185 kg (29 Stone)

Overall Dimension:

2130 x 900mm

Adjustable Platform Height:

67mm -640mm

Mattress Dimensions:

2000 x 900mm

Angle of Back Rest:

Up to 70 degrees

Angle of Knee Brake

Up to 30 degrees

Bed Weight:

81 kg

Key Information:

Manual Leg Rest Adjustment
Wooden Head and Foot Boards- Wrapped for decontamination
Low Voltage- 29V; 9V Battery Backup
Easy Cleaning- Complies to MHRA Guidelines
Comprehensive User Manual
CE Marked
Subject to the terms of Harvest Healthcare Guarantee and Warranties
the Ultra Low Bed is guaranteed for five years from date of purchase.

Supplier Contact Details:

Medguard Professional Healthcare

Tel: +353 (0) 1 835 2411

Email: sales@medguard.ie

Fax: 01 969 5009

Website: www.medguard.ie

Unit 6/7, Block 4, Ashbourne Business Park, Ashbourne, Co. Meath, A84 VK52