

# CORE PLUS



## High quality outdoor cabinet

- Protects against theft and vandalism due to pinlock and alarm
- Can withstand temperatures between -30°C and 50°C
- Dust and splash proof: IP55
- Able to house all types of AEDs
- Easily recognizable as AED cabinet due to appearance



### Dimensions

- 430 x 504 x 214 mm
- Weight: 4,25kg



### Features

- Audible alarm
- Internal lighting
- Two heating elements (24V)
- Digital keypad lock
- Four ventilators (24V)



### Delivery

- Worldwide shipping
- Shipment within 24 hrs