

M83

Suture Practice Trainers -Full set-

Various types of skins to improve the basic of suturing skills



FEATURES

- 1 | The artificial skin can be cut easily at any site to prepare a novel wound.
- 2 | The soft skins are made of durable materials.
- 3 | Simple and durable supporters keep the soft skins arched or maintain tension in the tubular models.
- 4 | All instruments for the practice session are included in the package.


SKILLS

- 4 skills for suture: abdominal skin suture, blood vessel suture, intestinal suture, limbs skin suture


KEY FEATURES




Abdomen



Blood vessel



Intestine



Limbs

DESCRIPTIONS

SPECIFICATIONS

SET INCLUDES

- | | |
|--|-----------------------------------|
| 1 set of simulated skins, ducts and supporters | 1 set of instruments |
| - abdomen set: 2 abdominal skins, 1 supporter | - needle-holders - suture needles |
| - blood vessels: 5 narrow tubes | - scissors - suture threads |
| - intestine set: 5 wide tubes, 1 supporter | - tweezers |
| - limbs set: 2 narrow skins, 1 supporter | 1 storage case |

MATERIALS

Soft resin
Latex free

RECOMMENDED DEVICES

REPLACEMENT PARTS

- | | |
|-----------|---------------------|
| 11219-030 | 1 abdominal skin |
| 11219-050 | 5 blood vessels |
| 11219-070 | 5 intestinal tracts |
| 11219-010 | 5 limbs skins |