


MW3

# Epidural Anesthesia Simulator

## Practice and perform epidural anesthesia with confidence

Epidural Anesthesia Simulator offers multiple training methods on both thoracic and lumbar areas, including lumbar punctures.




### FEATURES

- 1 | True-to-life resistance with the "pop" sensation
- 2 | Left and right lateral positions
- 3 | Simulates loss of resistance to saline or water
- 4 | The broad puncture pad allows paramedian access as well
- 5 | as median approach.

### SKILLS

- | Epidural anesthesia procedures at lumbar and thoracic area
- | Lumbar puncture



Landmark palpation



Epidural anesthesia



Lumbar puncture

### DESCRIPTIONS

#### SPECIFICATIONS

Size: W48 x D23 x H29 cm  
W18.9 x D9 x H11.4 inch

Weight: 5kg / 11lbs

#### SET INCLUDES

- 1 | life size unisex torso
- 5 | puncture pads
- 1 | supporting stand
- 1 | 50ml syringe
- 1 | drain pouch
- 1 | instruction manual

#### MATERIALS

Soft resin  
Latex free

#### REPLACEMENT PARTS

11382-010 2 replacement puncture pads

#### RECOMMENDED DEVICES

Epidural needle: 18G  
Spinal needle: 21G