

MW50

## Strap-on Venipuncture Simulator

### Cost-effective training unit for venipuncture

Strap-on Venipuncture Simulator is a simple, economical training kit for intravenous injection. This trainer is to be strapped onto the upper arm of an SP (Standardized Patient) or a life-size training manikin.



#### FEATURES

- 1 | Realistic tissue and vein wall resistance
- 2 | Vein tube can be filled and refilled easily using a soft plastic bottle that works as a manual pump.
- 3 | Three vein tubes are embedded: 2 normal and 1 thin.

#### SKILLS

- | Blood collection and intravenous injection

#### DESCRIPTIONS

##### SPECIFICATIONS

Size: W12 x D14 x H2.5 cm / W4.7 x D5.5 x H0.98 inch  
W6.5 x D11 x H2.0 cm / W2.6 x D4.3 x H0.79 inch (Pad)

##### SET INCLUDES

- 5 | puncture pads (vein tube embedded)
- 5 | supporting frames
- 5 | belts
- 5 | injection bottles
- 1 | instruction manual
- 1 | storage case

##### MATERIALS

Soft resin  
Latex free

##### REPLACEMENT PARTS

11431-010 5 replacement puncture pads

##### OPTIONAL PARTS

11430-010 Arm Vessels Anatomy Sleevelet  
11388-400 10 cotton swab red

##### RECOMMENDED DEVICES



-

M94

## Intradermal Injection Training Simulator

### Wheal can be shown by correct injection

This simulator provides training for intradermal injection to the inner forearm. Swelling of the skin will form when proper injection procedures are followed.



#### FEATURES

- 1 | Life-like touch and resistance of the injection site resembles the human arm
- 2 | Injection pads are easily changeable.
- 3 | No needle trace remains.

#### SKILLS

- | Intradermal Injection

#### DESCRIPTIONS

##### SPECIFICATIONS

Size: W9 x D20 x H8 cm Weight: 0.3 kg/ 0.7 lbs  
W3.5 x D7.9 x H3.1 inch

##### SET INCLUDES

- 5 | injection pads
- 1 | forearm model
- 1 | supporter stand
- 4 | injection pad cases
- 1 | instruction manual

##### MATERIALS

Soft resin  
Latex free

##### REPLACEMENT PARTS

11346-010 1 replacement puncture pad

##### RECOMMENDED DEVICES

Needle 29G

\*Needle not included.