

MW29




NEW

# Oral Care Simulator Advanced

## Oral care simulator which can mimic various patient conditions

True-to-life oral cavity and denture sets  
Actual water can be used to clean oral cavity



### FEATURES

- 1 | Viewable and accessible trachea and esophagus
- 2 | Interchangeable denture with tooth diseases
- 3 | Simple set-up and head positioning
- 4 | Includes washable plaque and oral thrush
- 5 | Suction training with water

### SKILLS

- | Oral cavity assessment
- | Care of the oral mucous membrane
- | Brushing and cleaning of the oral cavity
- | Identifying oral thrush
- | Care for humidity retention
- | Replacement and maintenance of the denture
- | Understanding of airway suction
- | Massage of oral cavity

### KEY FEATURES

#### Oral anatomy of denture with problems

- | spaced arch
- | food residue
- | anodontia
- | stomatitis
- | isolated tooth
- | oral thrush
- | crowding
- | gingival recession
- | root fracture
- | stump of a tooth



#### Patient conditions

- | Complete denture
- | Dental problems
- | Toothless
- | Intratracheal intubation
- | Tracheostomy
- | Tube feeding

#### Patient positions

- | Sitting position
- | Half sitting position
- | Lateral position

### DESCRIPTIONS

#### SPECIFICATIONS

Size: W21 x D37 x H26 cm / W8.3 x D14.6 x H10.2 inch Weight: 2.6 kg / 5.7 lbs

#### SET INCLUDES

1 head/neck manikin	1 simulated residue	1 storage case
1 base	1 mirror	1 instruction manual
1 normal denture	2 sponge-tipped swabs	
1 denture with tooth disease	1 toothbrush	
1 denture case	1 denture brush	

#### MATERIALS

Soft resin  
Latex free

#### REPLACEMENT PARTS

11408-020	normal denture
11408-030	denture with tooth disease
11408-010	simulated residue

\* Instruments are not included.