

MW43




"AbdoAbby" Abdominal Examination Simulator

True-to-life tactile representation of abdominal organs to cultivate the sense of touch and serve as a benchmark during palpation.

Product Supervision

Takahiro Amano, M.D.,
Senior Advisor for Academic Affairs and Professor
International University of Health and Welfare,
Honorary Director Sanno Medical Center, Japan

Kyoko Arai
RN, PHN, PhD, Professor,
School of Nursing, University of Shizuoka, Japan



FEATURES

- 1 | Supple, resilient, and delicate textures which allow both shallow and deep palpation
- 2 | Interchangeable organs
- 3 | Respiratory motion
- 4 | Abdominal auscultation

SKILLS

- | Visual Inspection
- | Auscultation
 - renal artery
 - abdominal aorta
 - iliac artery
 - bowel sounds
- | Palpation
 - liver, spleen, kidney, uterus and lymphoma
- | Percussion
 - liver, spleen and rib

KEY FEATURES

Bony landmarks


- | Pelvis
- | Ribs
- | Lower rib cage
- | Costal margin
- | Xiphisternum
- | Pubic crest
- | Anterior superior iliac spines

Liver

- | Normal
- | Chronic hepatitis
- | Early late cirrhosis
- | Cirrhosis

Inguinal area

- | Normal
- | Lymphoma



Spleen

- | Normal
- | Enlarged
- | Slightly enlarged

Kidney

- | Normal
- | Hydronephrosis
- | Polycystic kidney

Uterus

- | Normal
- | Myoma

Bruits and sounds

- Renal artery
- | Abdominal aorta
- | Iliac artery
- | Bowel sounds

DESCRIPTIONS

For product codes of the alternate skin tone versions of this model, please refer to www.kyotokagaku.com.

SPECIFICATIONS

Size: 40 x 80 x 20 cm / 15.7 x 31.5 x 7.9 inch
Weight: 14kg / 30.9lbs

MATERIALS

Soft resin
Latex free

SET INCLUDES

1	torso unit	2	inguinal lymph nodes (normal, lymphoma)
1	talcum powder	2	uterus (normal, fibroid)
1	set of the respiratory motion control unit	1	simulated stethoscope
2	kidneys (hydronephrosis, cysts)	1	controller
2	spleens (slightly enlarged, enlarged)	1	AC adapter
4	livers (precirrhosis, cirrhosis, chronic hepatitis, normal)	3	storage cases
1	talcum powder	1	instruction manual