




Nasopharyngeal Swab Collection Simulator

For training to healthcare professionals who conduct virus testing for the first time



FEATURES

- 1 | When swab reaches the right position in the nasal cavity, the swab will be pigmented in blue
- 2 | Ways to use tongue depressor with swabs can be practiced.
- 3 | The shape of nasal cavity between right and left are different.
- 4 | The product comes with cross sectional model.



SKILLS

Training to healthcare professionals who conduct virus testing for the first time.



DESCRIPTIONS

SPECIFICATIONS

Size: W42xD23xH38cm/W16.8x D9.2x H15.2in Weight: 2.15 kg / 4.73 lbs

MATERIALS

Soft resin, Hard resin
Latex free

SET INCLUDES

- 1 head manikin
- 1 cross sectional model
- 1 syringe (10ml)
- 1 disposable cup
- 10 cotton swab (blue)
- 6 pad (1 installed in the head manikin)
- 1 storage case
- 1 instruction manual

REPLACEMENT PARTS

- 11426-010 5 replacement pads
- 11402-010 10 cotton swab (blue)

