



# MIAMI J®

Known as the superior c-spine immobilizer through multiple independent studies, the Miami J's patented design meets the specialized needs of all those afflicted with cervical injuries. The Miami J offers features that enhance compliance while providing immobilization, inhibition of skin breakdown and patient comfort.



SPINE SOLUTIONS

## INDICATIONS FOR USE

- Post-trauma
- Stable cervical fractures
- Cervical spondylitis
- Cervical spondylosis
- RA and OA cervical spine
- Herniated cervical disc
- Post-operative
- Motor Neurone Disease

## FEATURES AND BENEFITS

- Sorbatex padding is antibacterial and clinically proven to inhibit microbial growth, which enhances comfort and helps assure patient compliance
- Bioengineered to minimize pressure points in key known areas: Chin, occiput, trapezius and clavicle
- Angulation buttons for occiput and sternal pad is removable when supine, lying and swallowing
- MRI, CT and X-ray lucent, MR safe

## REQUIRED MEASUREMENTS

- Sizing based on 6 phenotypes to ensure correct fit and optimum immobilization

Phenotype-driven sizing system makes size selection simple while drastically reducing margin of error in fitting a patient. Each collar has a bioengineered front with a specific back designed to work together for superior immobilization. For pediatric sizing, please see Miami Jr in this Catalog.

### Specialty Sizes (40% of population)

### Standard Sizes (60%)

Kyphotic (10%)	XSmall (5%)	Stout (15%)	Tall (10%)	Regular (30%)	Short (30%)
<ul style="list-style-type: none"> <li>• kyphotic</li> <li>• chin-on-chest</li> <li>• geriatric</li> <li>• osteoporitic and/or ankylosing spondylitic</li> </ul>	<ul style="list-style-type: none"> <li>• Very short, thin neck</li> <li>• Principally female</li> <li>• More prevalent in the Asian population</li> </ul>	<ul style="list-style-type: none"> <li>• Very large neck circumference</li> <li>• Obese</li> <li>• "No neck"</li> <li>• Massive shoulders</li> </ul>	<ul style="list-style-type: none"> <li>• Long, tall neck</li> <li>• "Swan neck"</li> <li>• Young Women</li> <li>• Adolescents</li> </ul>	<ul style="list-style-type: none"> <li>• Female</li> <li>• Standard size for mature women</li> <li>• Thin, mature male</li> </ul>	<ul style="list-style-type: none"> <li>• Standard size for mature men</li> <li>• Short-necked women</li> </ul>

\*For individuals with neck circumference > 25", the Bariatric collar is recommended.

# MIAMI J®

## ORDERING INFORMATION

### MIAMI J

Size	Collar set*	Collar	Collar Back	Collar Front
Super Short	MJR-200S	MJ-200S	MJ-400BK	MJ-200SFRT
Stout	MJR-200L	MJ-200L	MJ-200LBK	MJ-200LFRT
XS	MJR-250	MJ-250	MJ-250BK	MJ-250FRT
Short	MJR-300	MJ-300	MJ-300BK	MJ-300FRT
Regular	MJR-400	MJ-400	MJ-400BK	MJ-400FRT
Tall	MJR-500	MJ-500	MJ-500BK	MJ-500FRT
Bariatric	MJR-200B	MJ-200B		

\*Collar sets come with one additional set of replacement pads.

Miami J Front and Back Panels may be interchanged to accommodate patient anatomy, e.g. a Short Front may be used with a Tall Back. Choose collar size based on patient neck height, not patient body size. Super Short collars are packaged with Regular Back Panels

## ACCESSORIES

### REPLACEMENT PADS

Size	Part#
Super Short	MJP-100
Stout	MJP-200L
XS	MJP-250
Short	MJP-100
Regular	MJP-100
Tall	MJP-100
Bariatric	MJP-200B

### MIAMI J PATIENT PACK

Description	Part#
Miami J Patient Pack	MJ-PTPK2