



Jura Three Panel Immobiliser

A universal three panel knee immobiliser manufactured from a quality foam laminate with terry towel lining. This brace is designed to offer a lower cost alternative to hospital departments with high volume usage.

Can be used for an array of injuries and conditions such as:

- + Conservative treatment for knee injuries
- Maintain knee extension following surgery, trauma or soft tissue injuries



Code - 0°	Code - 20°	Splint Length
JTI/12	JTI/12/20D	30cm
JTI/14	JTI/14/20D	36cm
JTI/16	JTI/16/20D	41cm
JTI/18	JTI/18/20D	46cm
JTI/20	JTI/20/20D	51cm
JTI/22	JTI/22/20D	56cm
JTI/24	JTI/24/20D	60cm
JTI/26	JTI/26/20D	66cm
JTI/28	JTI/28/20D	71cm