



Choose the right stethoscope to help you perform at your best.

3M™ Littmann® Stethoscopes feature adjustable headsets with soft-sealing eartips. Select Littmann mechanical stethoscopes are designed with tunable technology to give you the ability to hear different frequency sounds by simply adjusting the pressure applied to the chestpiece. Get the acoustic performance, high quality and comfort you need with a Littmann stethoscope.



| | Lightweight II S.E. Stethoscope | Classic II Paediatric and Infant Stethoscopes | Classic III™ Monitoring Stethoscope | Cardiology IV™ Diagnostic Stethoscope | Master Cardiology™ Stethoscope | CORE Digital Stethoscope |
|--|--|---|---|---|---|--|
| | Ideal for taking blood pressure readings and making limited patient assessments. | Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another. | | Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation. Useful in challenging and critical care environments. | | High-performance auscultation. Hear more with amplification, noise cancellation and mobile connectivity. |
| Approximate Length | 71cm | 71cm | 69cm | 69cm | 69cm | 69cm |
| Approximate Weight | 118g | 105g, 95g | 150g | 177g | 185g | 232g |
| Warranty | 2 years | 3 years | 5 years | 7 years | 7 years | 2 years |
| Diaphragm | Tunable | Non-Tunable (Floating) | Tunable (2) | Tunable (2) | Tunable | Tunable (2) |
| Approximate Chestpiece Diameter | Double-sided 5.2 x 4.5cm 3.5 x 3.2cm | Double-sided Paediatric 3.4cm 2.7cm Infant 2.8cm 1.7cm | Double-sided Adult side 4.3cm Paediatric side 3.4cm | Double-sided Adult side 4.3cm Paediatric side 3.4cm | Single-sided 4.8cm | Double-sided Adult side 4.3cm Paediatric side 3.4cm |
| Features | Teardrop-shaped chestpiece Fits under blood pressure cuff | Specially designed for the smallest patients | Convertible open bell | Convertible open bell Dual-lumen tubing | Dual-lumen tubing Best acoustic response in Littmann mechanical line | Up to 40x amplification* Switch from analog listening to digital Connects to Eko software on a smart device** to capture, store and share data |



The power to learn, in the palm of your hand.

3M™ Littmann® Learning App – Download free version in the App Store® or Google Play™ Store

*At peak frequency, vs. analogue mode
**Not included