

Universal Mounting Brackets

The universal accessory pack contains all the necessary parts to allow the Sharpak 18, 35, 60, 120 and 220 PLUS containers to be safely mounted either on a wall or fixed to a trolley or equipment stand.




Patented and trademarked by Nolato.

PRODUCTS

PRODUCTS	CODE	QTY
Wall & trolley bracket for 18 PLUS	SM06500	1
Wall & trolley bracket for 35 - 220 PLUS	SM05000	1



WALL, STAND OR TROLLEY MOUNT FOR 18 - 220 PLUS



STAND OR TROLLEY MOUNT FOR 18 - 220 PLUS



WALL MOUNT FOR 18 - 220 PLUS





Universal Wall bracket fixing height

Where possible it is recommend that sharps containers should be located in a safe place fixed to wall brackets.

In order to ensure that the container aperture is at the best working height, whether measured from the floor or from the work surface, the bracket can be correctly positioned by using the simple formula shown below.

Wall and trolley brackets can be supplied free of charge. Please contact us for more information.

