

TERi™ Androgynous Geriatric Trainer



Product Name	TERi™ Androgynous Geriatric Trainer	
Addresses Training Needs/Curricula	<p>General capability carefully selected to support practice on:</p> <ul style="list-style-type: none"><li>- Essentials of Patient Care</li><li>- Movement Assistance and Fall Prevention</li><li>- Grooming and Daily Living Assistance</li><li>- Basic and Advanced Nursing Skills</li><li>- Cardiovascular and Respiratory Emergencies</li></ul> <p>Suitable for long-term facilities and emergency room training.</p>	
Configuration	<ul style="list-style-type: none"><li>• Light or medium skin tones.</li><li>• Configurable to fit unique training needs and budgets.</li><li>• A range of accessories available that enhances the baseline version maximizing versatility, training benefit and return on investment.</li></ul>	
Accessories	<p>IV Arm (Light) - LF04303 IV Arm (Medium) - LF04303M</p> <p>BP Arm (Light) - LF04304 BP Arm (Medium) - LF04304M</p>	
Consumables	<ul style="list-style-type: none"><li>• IV Arm:<ul style="list-style-type: none"><li>- Skins</li><li>- Veins Kit</li><li>- IM Pads</li><li>- Fluid Supply Bag</li></ul></li><li>• Urinary Catheterization:<ul style="list-style-type: none"><li>- Male Genitalia Foreskins</li><li>- Foley Catheter</li><li>- Urinary Catheterization Bladders</li></ul></li></ul>	<ul style="list-style-type: none"><li>• IM Pads for Thighs and Buttocks</li><li>• Lung Bags</li><li>• Stomach Bags</li><li>• Stoma Bags</li><li>• Gastrostomy Bags</li><li>• Simulated Blood</li><li>• Lubricant</li></ul>
Product Code/SKU	<p>Patient Care Trainer (Light) - LF04300 Patient Care Trainer (Medium) - LF04300M</p> <p>Patient Skills Trainer (Light) - LF04301 Patient Skills Trainer (Medium) - LF04301M</p> <p>Auscultation Patient Skills Trainer (Light) - LF04299 Auscultation Patient Skills Trainer (Medium) - LF04299M</p> <p>Patient Simulator (Light) - LF04302 Patient Simulator (Medium) - LF04302M</p>	
Dimensions/Weight	<p>Dimensions: 5 feet Ship Weight: 65 lbs or 29.48 kg</p>	



Nasco Healthcare Inc.  
16 Simulaids Drive, Saugerties, NY, 12477 USA  
Phone: 1-833-NASCOHC (627-2642) Email: info@nascohealthcare.com  
www.nascohealthcare.com



TERi™ Androgynous Geriatric Trainer





Be **READY** when called into action.

**TERi™ Androgynous Geriatric Trainer**

A full-body elderly patient care trainer configurable in male or female variations for medical simulation. Evolving from our GERi™/KERi™ manikins, TERi™ has a more realistic look and feel to represent human anatomy, with weight distributed to represent a real patient for lifting and carrying.

Each version is affordable and durable with upgradable items.

- **Patient Care Trainer<sup>1</sup>** – General patient care, daily living assistance simulation platform.
- **Patient Skills Trainer<sup>2</sup>** – Physical skills practice simulation platform with injection/IV, catheterization, cardiovascular and respiratory support.
- **Auscultation Patient Skills Trainer<sup>4</sup>** – All the features of Patient Skills Trainer with added auscultation. Lung sounds at five anterior, ten posterior, and two midaxillary locations. Six anterior heart sites.
- **Patient Simulator<sup>3</sup>** – Comprehensive and complete simulation platform including patient monitoring capability.



The next generation of elderly care simulation is here.



**PROM: Passive Range of Motion<sup>1-4</sup>**  
Seating, Laying, Finger/Toe Manipulation with No Effort from Patient

**Cardiovascular Support<sup>2-4</sup>**  
Carotid Pulses, CPR with Quality Performance Metrics & Auscultation



**Intravenous, Injection, Subcutaneous Injection Sites<sup>2-4</sup>**, optional on 1  
Multiple Locations Equipped with Standard Arm IM Injection Pads



**Nasogastric Tube Procedures<sup>2-4</sup>**  
Placement, Lavage, Gavage, Feeding and Suctioning




**Enhanced Aesthetics<sup>1-4</sup>**  
Soft Silicone Skin with Noticeable Normal and Cancerous Skin Moles, Stage 1 Sacral Ulcer and Reddened Skin Folds

**Respiratory Support<sup>2-4</sup>**  
Manual Ventilation with Chest Rise, Oxygen Administration, Tracheostomy Care, Airway/Trach Suctioning, Airway Intubation, Ventilation and Auscultation






**Grooming & Daily Living Assistance<sup>1-4</sup>**  
Bathing, Clothing Changes, Denture Placement/Removal, Hearing Aid Placement/Removal, Eye Irrigation & Optic Drops, Ostomy & Urinary Catheter Care, Hair Care



**Indwelling Male & Female Urinary Catheterization<sup>2-4</sup>**  
Includes Support for Urethral and Suprapubic Catheters

**Advanced Skills<sup>2-4</sup>**  
Urinary Catheterization, Enema Administration, Pelvic and Prostate Examination



Available in light or medium skin tones.  

This manikin is intended to be used as an instructional aide for teaching purposes. Actual product may vary slightly from photo. Nasco Healthcare reserves the right to change product color, materials, supplies, or function as needed.