

Data Sheet

THERMARMOUR Baby Wrap

Product Code	OB202-55-57---TAEM55
GTIN/UDI	05060885279682
Colour(s)	Emergency Blue Outer with White Inner
Size	55cm x 57cm (can be cut to size if required)
Material	Non-conductive bilaminate foil & needle punched polyester fleece
Care	Patient Specific, cannot be laundered
Buy Online	www.interweavetextiles.com/products/patient-warming/

Designed with the goal of maintaining and restoring normothermic temperatures in emergencies, THERMARMOUR is suitable for patients at risk of hypothermia. Other scenarios include, but are not limited to, pre-warming for surgeries and temperature management of trauma casualties.

Key Features	<ul style="list-style-type: none"> ❖ Latex Free, Recyclable & X-ray Compatible. ❖ Easily transportable standalone device. ❖ No waiting to use once opened; a bonus for first responders. ❖ A cost and time efficient passive warming method.
---------------------	--

Usage

Open packaging and remove the THERMARMOUR Baby Wrap. Unfold the wrap fully, placing the infant in the middle as shown. Wrap and secure with the foil side outwards.



Packing

Manufactured in our UK facility, THERMAMOUR is a readily available, sustainable solution.


Per Carton	40, Individually Bagged
Carton Size	40 x 32 x 30cm
Pallet Quantity	2160

Data Sheet

Storing and Handling

Keep dry and store out of direct sunlight until needed. THERMARMOUR products are for single patient use and should be discarded once soiled. If there is any damage identified upon receipt, please contact sales@thermarmour.com before use.

Technical Drawing



Risks and Precautions

Following the trial and testing period, no known residual risks or side effects have been identified. Please ensure all components including packaging are recycled or disposed of safely, as misuse may cause suffocation.