

| Name                                                                                  | Value                                                               |
|---------------------------------------------------------------------------------------|---------------------------------------------------------------------|
| Pharmacy refrigerator                                                                 | MD14DD-P301QX23AXXXZ1SN4B                                           |
| EAN code                                                                              | 4770104014048                                                       |
| Gross capacity, l                                                                     | 130                                                                 |
| Net capacity, l                                                                       | 100                                                                 |
| Factory temperature setting, °C                                                       | +5                                                                  |
| Fluctuation, °C                                                                       | 0,5                                                                 |
| Gradient, °C                                                                          | 2,9                                                                 |
| Increase in temperature from +5 °C to +10 °C (empty, +25 °C TU)                       | 40 min.                                                             |
| Decrease in temperature from +25 °C to +5 °C (empty, +25 °C TU)                       | 40 min.                                                             |
| Recovery time after 1 min of door opening (empty, +25 °C TU, according to EN 60068-3) | 8 min.                                                              |
| Positive temperature deviation from the setpoint                                      | 0.5                                                                 |
| Negative temperature deviation from the setpoint                                      | -1.4                                                                |
| Heat output, kJ/h                                                                     | 135                                                                 |
| Ambient temperature range, °C                                                         | from +10 to +35                                                     |
| Airborne acoustical noise emissions, dB(A) re1 pW                                     | 42                                                                  |
| Refrigerant                                                                           | R600a                                                               |
| Daily consumption, kWh/24h                                                            | 0,751                                                               |
| Annual energy consumption, kWh/year                                                   | 274                                                                 |
| Rated voltage, V                                                                      | 220 - 240                                                           |
| Rated frequency, Hz                                                                   | 50                                                                  |
| Access port                                                                           | Yes                                                                 |
| Calibration block                                                                     | Yes                                                                 |
| Volt-free alarm contact                                                               | Yes                                                                 |
| Self closing door                                                                     | Yes                                                                 |
| Cooling system                                                                        | Fan forced cooling                                                  |
| Defrost                                                                               | Automatic                                                           |
| Interior light                                                                        | LED ceiling lighting                                                |
| Malfunction warning signal                                                            | Optical and acoustical                                              |
| Power failure alarm                                                                   | Immediately upon power failure for min. 12h                         |
| Connection to the external alarm system                                               | Yes                                                                 |
| Alarm indicator on screen                                                             | Visible until removed manually                                      |
| Data logger                                                                           | Mobile APPs connected via Bluetooth                                 |
| Protection against low temperature                                                    | Yes                                                                 |
| Number of shelves                                                                     | 2 adjustable + 2 bottom, wire shelves coated by white polymer paint |
| Shelf size (W x D), mm                                                                | 420 x 324, bottoms 465 x 263 and 445 x 180                          |
| Maximum load per shelf, kg                                                            | 20                                                                  |

| Name                                     | Value                                             |
|------------------------------------------|---------------------------------------------------|
| Shelf material                           | Wire shelves coated by white polymer paint        |
| Side wall material                       | Steel, white painted                              |
| Inner cabinet material                   | HIPS polystyrene                                  |
| Type of door                             | Double panel selective glass in aluminium profile |
| Type of lock                             | Mechanical                                        |
| Foot                                     | 2 adjustable foot at the front, 2 back            |
| External dimensions (H x W x D), mm      | 870 x 560 x 620                                   |
| Internal dimensions (H x W x D), mm      | 725 x 465 x 392                                   |
| Package dimensions (H x W x D), mm       | 950 x 570 x 640                                   |
| Weight net, kg                           | 38                                                |
| Weight gross, kg                         | 41                                                |
| Pharmacy refrigerator packed on a pallet | MD14DD-P301QX23AXXXZ1SN4A                         |
| EAN code                                 | 4770104014048                                     |
| Dimensions (H x W x D), mm               | 1080 x 570 x 710                                  |
| Weight gross, kg                         | 49,5                                              |