

Name	Value
Pharmacy refrigerator	MD35DD-P301DDQMAXXX1SN4B
EAN code	4770104035050
Gross capacity, l	350
Net capacity, l	240
Factory temperature setting, °C	+5
Fluctuation, °C	0,5
Gradient, °C	2,9
Increase in temperature from +5 °C to +10 °C (empty, +25 °C TU)	40 min.
Decrease in temperature from +25 °C to +5 °C (empty, +25 °C TU)	40 min.
Recovery time after 1 min of door opening (empty, +25 °C TU, according to EN 60068-3)	8 min.
Positive temperature deviation from the setpoint	1.3
Negative temperature deviation from the setpoint	-1.1
Heat output, kJ/h	266
Ambient temperature range, °C	from +10 to +35
Airborne acoustical noise emissions, dB(A) re1 pW	50
Refrigerant	R600a
Daily consumption, kWh/24h	1.776
Annual energy consumption, kWh/year	648
Rated voltage, V	220 - 240
Rated frequency, Hz	50
Self closing door	No
Access port	Yes
Calibration block	Yes
Volt-free alarm contact	Yes
Cooling system	Fan forced cooling
Defrost	Automatic
Volt-free alarm contact	Yes
Interior light	LED ceiling lighting
Malfunction warning signal	Optical and acoustical
Power failure alarm	Immediately upon power failure for min. 24h
Alarm indicator on screen	Visible until removed manually
Data logger	Mobile APPs connected via Bluetooth
Protection against low temperature	Yes
Number of shelves	5 adjustable + 1 bottom
Shelf size (W x D), mm	523 x 351, bottom 445 x 180
Maximum load per shelf, kg	40

Name	Value
Shelf material	Wire shelves coated by white polymer paint
Side wall material	Steel, white painted
Inner cabinet material	HIPS polystyrene
Type of doors	2 double panel selective glass in aluminium profile
Type of lock	Mechanical
Adjustable foot at the front	2
Wheels behind	2
External dimensions (H x W x D), mm	1810 x 600 x 620
Internal dimensions (H x W x D), mm	1518 x 502 x 453
Package dimensions (H x W x D), mm	1870 x 620 x 630
Weight net, kg	77
Weight gross, kg	82
Pharmacy refrigerator packed on a pallet	MD35DD-P301DDQMAXXXX1SN4A
EAN code	4770104035067