

SKID OK-B

Evacuation chair – Silver frame



Skid-Ok is the new ultra compact Spencer evacuation chair with a new aluminum frame and tubular trellis swing fork which characterises the design. Technically advanced, Spencer Skid-Ok represents the highest standard of quality.

Manufacturing materials and components are all innovative and contribute to make the chair lightweight and easy to manoeuvre which means it is ready to offer great performance in any scenario.

The unquestionable comfort of the device is underlined by reclining armrests with rubber finish.

The innovative self-braking slide system offers the operator total control during the stair descent.

Manufactured in welded aluminium with polished finish.

Adjustable seat strengthened with Spentex®.

Ø100 x 32 Front Pivoting wheels with break and parking break Ø 150 x 30 Rear wheels

Skid -OK is supplied with 2 x 50 mm polypropylene belts.

Technical details

Width	520 mm
Length	900 mm
Height with big handles open	1540 mm
Height with big handles closed	1040 mm
Width closed	175 mm
Weight	10 kg
Max loading weight	150 kg
Surface treatment.	Polishing
Front wheels n°2	Ø 100 x 32 mm
Posterior wheels n°2	Ø 150 x 30 mm
Canvas	Spentex ®
Canvas color	Black/Red
Frame	Aluminum
Frame color	Chromed
Belts	50 mm polypropylene
Latex free	
SK20002E	SKID OK-B

