

# MS LED MS LED PLUS

Minor surgery luminaire and advanced exploration with LED light source and articulated arms of great precision allow for regulation of the luminaire in any position.

The quality of the materials used in its manufacture makes it especially suitable for intensive use in surgery and examinations.

## Exceptional quality of LED light

The powerful LED light source offers exceptional light quality and a lifespan of more than 50,000 hours.

The high performance of the colour Ra 97% depicts the exceptional quality of white light for an excellent distinction of fabrics and reproduction of all visible colours, naturally and precisely.

## Versions

<b>MS LED</b>	25.000 lux at 50 cm On/Off version
<b>MS LED PLUS</b>	45.000 lux at 50 cm Adjustable light intensity version

## Ergonomic design

The circular shape of the head and its handle in the lower part facilitate its exact positioning, fulfilling the ergonomic requirements for any use.

Designed with clean and soft shapes which facilitate hygiene and cleanliness in hospital environments.

## Mounting options

Standard version includes AH table clamp.

As accessories, we offer a wide selection of supports such as wall clamps, wall extension arm and trolley stand.



●..... MS LED Plus version.  
Dimmable, adjustable light intensity.

MS LED version



## MS LED

Technical Details	(Cod. M183500)
Light source	LED 8W
Lux a 50 centimeters	25.000 lux
Adjustable light intensity	No
Colour temperature	4.500 °K
Colour rendering	> 97%
Lamp lifetime	50.000 hours
Arm length	90 cm
Colour	White

## MS LED PLUS

Technical Details	(Cod. M187001)
Light source	LED 12W
Lux a 50 centimeters	45.000 lux
Adjustable light intensity	Yes
Colour temperature	4.500 °K
Colour rendering	> 97%
Lamp lifetime	50.000 hours
Arm length	90 cm
Colour	White



**Wall support**  
Cod. B0101800



**AH table clamp**  
Cod. B010190



**Rail Plus bracket**  
Cod. 82026



**Extension wall arm**  
Cod. B0101S0

**Trolley stand**  
8,8 kg Cod. 09651

