

HEALTHCARE PEDAL BIN RANGE

Large Healthcare Handsfree Pedal Bin



PRODUCT CODE
2080

DIMENSIONS

Height: 81cm
Width: 38cm
Depth: 40cm

Volume: 70L
Bag: 26' X 44'



Anti-Bacterial Paint Standard Finish



Certified I.I.R.S.
Ref. no. S11210/2

Fire Resistant to
BS476 Part 7 - Class 0

Silent Closing

Slow Close Unit
available to prevent
Lid Slamming

Special Coating

Double Pass Paint
System available to
Help Prevent Rusting

Codes Breakdown

Start with 4 Digit Code
Add

SC for **Slow Close**

W for **White**, **Y** for **Yellow** etc
P before Code for **Primer**

Medium Healthcare Handsfree Pedal Bins



PRODUCT CODE
2082

DIMENSIONS

Height: 72cm
Width: 33cm
Depth: 32cm

Volume: 50L
Bag: 20' X 37'



PRODUCT CODE
2083

DIMENSIONS

Height: 58cm
Width: 39cm
Depth: 40cm

Volume: 50L
Bag: 26' X 44'

Features

- Fire Retardant
Cert Ref. SII 210/2
- Anti-Bacterial Finish
- Medical Grade Non
Porous Shock Cord
- Stainless Steel Pedal
- Heavy Duty Lock
- Screening or Labels
- Castors available
- Rubber Wall Buffers
- Rubber Floor Buffers

Small Healthcare Handsfree Pedal Bins



PRODUCT CODE
2085

DIMENSIONS

Height: 57cm
Width: 21cm
Depth: 27cm

Volume: 20L
Bag: 15' X 24'



PRODUCT CODE
2020

DIMENSIONS

Height: 57cm
Width: 21cm
Depth: 27cm

Volume: 20L
Bag: 15' X 24'

**ALL COLOURS
AVAILABLE**



Range of Bags Available