

OLYMPIC 180 HOIST



Product Code:	HLA034
Description:	<ul style="list-style-type: none">• The Olympic 180 Hoist benefits from the rigorous ergonomic design procedures and quality checks.• It has been constructed from high quality powder coated steel that will provide years of reliable service.• The hoist boasts an impressive 180kg lifting capacity, powered leg opening and market leading Linak® actuators and motors.
Safe Working Load:	180kg (28 Stone)
Lifting Range:	1755mm
Product Length:	1265mm
Product Height:	1335mm/2010mm (highest position/lowest position)
Leg Inside Width:	490mm/1075mm (closed/open)
Castor Size:	75mm/100mm (Front/Rear)
Carrier Bar Type	3 Hooked safety attachments for loop slings
Key Features:	<ul style="list-style-type: none">• Cost effective solution without compromising on quality• Robust powder coated steel frame• Linak control system• Electric leg opening as standard• Red emergency lowering screw• Easy to position with 4 large swivel castors and ergonomic handle• Tool free assembly• Single person installation• Can lift a patient from the floor