



**Code:** 34582  
**Category:** Oxygen concentrators  
**Unit of sale:** 1 pc.  
**Minimum order:** 1  
**Type:** Medical device  
**Class:** II A  
**NSIS:** 2449622  
**CND:** Z12159004  
**EAN13:** 6945040108245

**Description:** Small and portable oxygen concentrator composed of air compressor, absorption tower of molecular sieve, sensor and flowrate controller.  
The machine adopts PSA (pressure swing absorbers) method to directly concentrate medical oxygen with high concentrator from the air. Lightweight, easy to use and safe. Stable performance and low noise level.  
Wide application in clinic, hospital, oxygen bar, gym center...

Supplied with humidifier bottle, power cord, nasal cannula and user manual (GB, FR\*, IT, ES\*, DE\*) \* only on [www.gimaitaly.com](http://www.gimaitaly.com)

Size: 508x260x h 530 mm.

Weight: 21 kg

**Technical Specifications:** Max. flowrate: 5l/min  
Flow range: 0.5~5l/min  
Oxygen concentration: 93% ±3%  
Output pressure: 20kPa~50kPa  
Input power: 450VA  
Noise level: =55dB(A)  
Operating voltage: AC220V ± 10%, 50Hz ±1 Hz  
Type of protection: Class II, Type BF, IP21