

OXYGEN CONCENTRATOR 10 L

Code: 34572
Category: Oxygen concentrators
Unit of sale: 1 pc.
Minimum order: 1
Type: Medical device
Class: II B
NSIS: 2582202
CND: Z12159004
EAN13: 6933257924451



Description: SILENT OXYGEN CONCENTRATOR 10 I

Medical grade, silent, small oxygen concentrators for home use.
Provide high oxygen concentration through nasal cannula.
Work continuously.

Supplied with filter, humidifier bottle, nasal cannula and manual GB, IT.
Available on gimaitaly.com: ES, SE, GR, TR, Arabic.

Technical Specifications: Flow range: 1 ~ 10 l/min
Oxygen concentration: 87 - 96 %
Output pressure: 150 kPa
Input power: 700 VA
Noise level (front - whole unit): 49 dBA
Size - weight: 41x37.5x66.8 cm - 24 kg
Packing size: 51x44x75 cm
Operating voltage: 230 V - 50 Hz
Optimal oxygen concentration: after 15 min.
Classification: Type BF, Class II equipment