

CRASH / EMERGENCY CARTS

MX Emergency / Crash Carts were designed with input from emergency room professionals to deliver the ultimate in convenience and safety. These carts define a new standard of organization and function with specialised accessories for every emergency response. This series offers stable, smooth mobility and enhanced convenience with changeable accessories, as well as a breakaway lock for that little bit extra.



SPECIFICATIONS

	MX28EMG		MX31EMG		MX34EMG		MX37EMG	
	mm	inch	mm	inch	mm	inch	mm	inch
W	690	27.2	690	27.2	690	27.2	690	27.2
D	518	20.4	518	20.4	518	20.4	518	20.4
H	848	33.4	927	36.5	1006	39.6	1085	42.7
H (w/o 5" Casters)	714	28.1	793	31.2	872	34.3	951	37.4



2 Defibrillator shelf is able to turn 360° and lock tight during transportation.

Breakaway lock 2 options



First drawer only



Locking bar for all drawers

STANDARD ACCESSORIES

- 1 IV pole
- 2 Defibrillator shelf
- 3 CPR board
- 4 Oxygen tank holder
- 5 Sharps container set



CPR board
Macromolecular polymer CPR board is strong and easy to clean.



Oxygen tank holder
Adjustable ties securely fasten the oxygen tank.



Sharps container set