


Product Specifications Sheet

Specification Sheet No	SSN-FT1580-001			
Product Name	Non-Woven Swabs			
Product Description	Non-Woven Swabs			
Range	Premier			
Manufacturer	Henan Zhende Medical Products Co Ltd			
Material	50% polyester 50% viscose			
Sterilisation	ETO			
Issue Date	12 th September 2024	Version No	R03	

Intended use	Single use non-woven swabs are intended to be used to absorb fluids for medical purposes. Used to cover/clean wounds or abrasions and absorb their exudates or to apply topical medications.
---------------------	--

Packaging						
Product code	Colour	Dimension (cm)	Ply	Items Per Pack	Qty- Inner	Qty – Shipper Carton
1850A	White	5 x 5 (±0.5)	4	5	48 bags of 30 pouches	1440
1860A	White	7.5 x 7.5 (±0.5)	4	5	15 bags of 40 pouches	600
1861A	White	7.5 x 7.5 (±0.5)	4	2	15 bags of 100 pouches	1500
1870A	White	10 x 10 (±1)	4	5	10 bags of 40 pouches	400
1871A	White	10 x 10 (±1)	4	2	10 bags of 100 pouches	1000
1884A	White	10 x 10 (±1)	4	5	10 bags of 30 pouches	300

Quality Certification	CE certificate G21 001594 0003 Rev.00		
Device Classification	Class Is	Rule	4
Conformity Route	Medical Device Regulation (MDR) EU2017/745		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 Years- if stored as recommended above		
Contraindications/Warnings	N/A		
Markings	Premier Logo, Size, Qty, Ref, Lot, Expiry, CE marking, Latex Free, Sterilisation Method, Single Use		
Additional information	These swabs have high levels of absorption, fast wicking and are extremely low linting. They also have low adhesion which is important when treating open wounds.		
<i>Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.</i>			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	16/06/2023	Update to Mediq UK Brand	Jack Wong
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong
R03	12/09/2024	Update to change from MDD to MDR	Elizabeth Ogunkolade

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq UK Ltd.