

M-420 Digital Portable Stand On Scale



- Class III and MDD Approved
- Accurate to 50g
- Lightweight and portable
- Ideal for community health visitors
- Suitable for GPs, hospitals and the home
- Optional carry case available

Marsden M-420:

Digital Stand On Scale

- Class III and MDD Approved
- Accurate to 50g
- Perfect for use in the community, surgery or hospital
- Powered by 6 x AA batteries
- 100 hours of continuous use
- Auto power-off function helps to save battery power
- Lightweight and highly portable
- Optional carry case available to aid portability
- Easy to read display
- Hold and Tare functions to enable accurate and easy weighing
- Made from reinforced aluminium structure
- Covered by the Marsden 8 year warranty



Simple weighing

The M-420 is a heavy duty stand on scale, made from reinforced aluminium, which is ideal for weighing adults as it is accurate to a minimum of 50g.

The scale features from a 'tap-on' function, which means that the user does not need to bend down to operate the scale, and you can simply tap the device with your foot and it will power on.

Ideal for the community

We understand that in this day and age, patients are weighed in their own home, just as frequent as they are within a medical environment.

This scale is heavy duty and portable, making it the perfect additional to the arsenal of any community health visitor. Plus, it is easy to operate by medical professionals or people in their own home.

Easy to use

With three simple buttons and a large easy to use display, the M-420 could not be simpler to use and this is why it is one of our best-selling products.

The On/Off button ensures that you have full control over the device. Hold enables you to gain weights even for the most active babies and the tare keys allows you deduct the weight of the sling or any clothing.

Capacity:	220kg	Power Supply:	Mains/6 x AA Batteries
Graduations:	50g<150kg>100g	Battery Life:	100 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Class III and MDD Approved
Dimensions:	320mm x 310mm x 60mm	Optional Extras:	Bluetooth, Carry Case
Weight of Scale:	4kg	Warranty Period:	8 Years