

PRIME COMFORT ACTIVE MATTRESS



Product Code:

HHMPRIMA

Description:

The Prime Comfort Active mattress combines our Prime Comfort Plus with an active air cell component, giving the following outstanding features:

- Firm 'U' shaped high density polyether foam base for roll prevention and to aid transfers
- Castellated soft and resilient central foam insert
- 10 alternating cells powered by the low noise level pump
- Fully welded, antibacterial, water resistant, four way stretch cover
- Durable PVC base
- Step up: step down with/ without pump

Combining benefits of both static foam and alternating air cells, the Prime Comfort Active is suited to patients at very high risk of developing pressure ulcers.

Dimensions:

2000 x 900 x 140mm

Pump:

Compatible with Harvest 1 and Harvest 250 Pumps

Pump Dimensions:

280 x 140 x 100mm

Mattress Weight:

14kg

Max. User Weight:

39 Stone / 247kg

Pressure Relief:

High (without pump)

Very High (with pump)

Cover:

Welded, four way-stretch, vapour permeable, PU coated cover in plum, zipped to a durable PVC coated base section. Silver-based antibacterial properties in the PU coating are wash proof and environmentally friendly. Cover does not contain latex.

Cleaning Guidance:

Cover should be wiped down with hot water and mild detergent.

If heavily soiled, the cover may be machine washed up to 65°C and tumble dried. For disinfection use sodium hypochlorite solution max 1000ppm. Use no bleach, biological or phenolic-based cleaners. Dry completely before storing.

Warranty

In terms of Harvest Healthcare Guarantee and Warranties the Prime Comfort Active mattress is guaranteed for three years from date of purchase.

Exclusions:

1. Removable covers are guaranteed for one year
2. Damage to fabric covers through incorrect use and penetration by sharp instruments is excluded
3. Components of pumps which are recommended for annual service, are guaranteed for one year.

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Harvest Healthcare Ltd supplies a range of products such as foam mattresses, active air mattresses, fall protection mats and bed rail protection bumpers. These products are usually covered in vapour permeable PU fabric that could cause risk of suffocation in the event of prolonged facial contact.

(Mattress risk categories are to be used as a guide only)

Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 3: May 2021