

SCLEROTHERAPY/FILLER LUER NEEDLES 32G 0,23x12 - pink



Product Code	23684
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2648226
CND	A019006
EAN/UPC	8023279236842

Description

SCLEROTHERAPY / FILLER LUER NEEDLES - STERILE

- Gauge 32G (Ø 0.23 mm)
- Length: 12 mm
- Colour: pink
- Box of 100 pcs.
- Sterile

Meso-relle disposable hypodermic needles are multiple sharpened to facilitate penetration. Stainless steel AISI 304 cannula and polypropilene connection for syringe.

Sterilized with ethylene oxide.
Sold in individual blister made of medical paper.

Made in Japan.