




WOBURN COMMUNITY 1200 PROFILING BED



Product Code:

HLBC791.12

Description:

The Woburn Community 1200 is Harvest's first wide profiling bed, built to the same high quality as other beds in the range. Although not a bariatric bed, it has a SWL of 30 Stone. The generous width allows patients to move freely. The lockable handset allows for quick adjustments of height, back rest and knee break and the attractive vinyl-clad head and foot boards facilitate infection control.

Safe Working Load:

193kg (30 Stone)

Max. Patient Weight:

158kg (25 Stone)

Dimensions:

2120 x 1320mm

Mattress Base:

2000 x 1200mm

Max. Mattress Depth:

205mm

Adjustable Height:

385 - 815mm

Angle of Back Rest:

Up to 70 degrees

Angle of Knee Break:

Up to 30 degrees

Bed Weight:

91.5kg (without headboards)

Key Information:

4 Section Profiling Bed
Electronically Powered Back Rest & Thigh Rest
Manual Leg Rest Adjustment
Vinyl Coated Head and Foot Boards
Integral Side Rails
Trendelenburg Function as standard
Lockable Handset
Low Voltage – 29v; 9v Battery Backup
Comprehensive Manual; CE marked