


## WOBURN ULTRA-LOW 800 PROFILING BED



**Product Code:**

HLB800.01AR

**Description:**

A fully functioning four- section Ultra-Low profiling bed. Platform can be lowered to 90mm off the floor for patient safety and raised to 800mm to make it easier for the carer when providing care for the patient. Electronically setting of height, back rest or knee brake is simple, using the lockable handset. The Woburn Ultra-Low 800 has a Trendelenburg function as standard. Optional bed levels, lighting and accessories available.

**Max Patient Weight:**

175kg (27 Stone)

**Safe Working Load:**

185kg (29 Stone)

**Overall Dimension:**

2399 x 900mm

**Adjustable Platform Height:**

90mm -800mm

**Mattress Dimensions:**

2000 x 900mm

**Angle of Back Rest:**

Up to 70 degrees

**Angle of Knee Brake**

Up to 30 degrees

**Bed Weight:**

125 kg

**Key Information:**

- Low voltage
- Lockable handset
- Emergency battery backup
- Wooden head and foot boards
- In-built auto-regression
- Trendelenburg as standard
- Electric backrest and kneebreak with mechanical lower leg support
- Easy Cleaning - complies to MHRA guidelines
- Compatible with environmental control system
- Comes with a 5-year warranty