

WOBURN LOW PROFILING BED



Product Code:

HLB797.03

Description:

The dependable 4 section profiling bed now features high quality control system, a reduced minimum platform height and a higher load capacity. By lowering the platform to only 220mm from floor level, the risk of injury for those with a tendency to roll out of bed is significantly reduced. Comfort setting of height, back rest and knee break angles is easy when using the lockable handset.

Safe Working Load:

220kg (34 Stone)

Max. Patient Weight:

185kg (29 Stone)

Mattress Base Size:

2000 x 900mm

Max. Mattress Depth:

205mm

Adjustable Height:

220 - 650mm

Angle of Back Rest:

Up to 70 degrees

Angle of Knee Break:

Up to 30 degrees

Bed Weight:

89.5kg (without side rails)

Key Information:

4 Section Profiling Bed
Electronically Powered Back Rest & Thigh Rest
Manual Leg Rest Adjustment
Wraparound Wooden Head and Foot Boards
Side Rails - Optional
Trendelenburg Function and Auto-Regression - Optional
In-Built Extension - Optional
Lockable Handset
Low Voltage – 29v; 9v Battery Backup
Comprehensive Manual
CE marked