



Product Code	28073
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	1757691
CND	Z12159002
EAN/UPC	8033717011436
PARAF	977812181

Description

MYNEB NEBULIZER - piston

Compact, fast and silent nebulizer.
Suitable for the treatment of the respiratory tract diseases, bronchial and asthmatic diseases, allergies.
Ideal for nebulizing any medicine, even oily ones such as cortisones essentials oils, and active ingredients with solid elements (suspensions).

Supplied with carrying bag, fasterjet bulb, adult and pediatric mask, air tube, mouthpiece and nosepiece.
Multilingual box and manual: GB, FR, IT, ES, DE.
Class IIa.
Made in Italy.

Technical Specifications

AC adapter:
- input 100-240 V - 50/60 Hz - 0,6 A
- output 12 V - DC 2A
Operating cycle: 30/30 min. ON/OFF
Noise level: 51 dB
Nebulization rate: 0.45 ml/min with NaCl 0.9%
Air Flow: 12 l/min
Max pressure: 165 kPa/1.65 bar
Particle size: 75% < 5 µm (NaCl 0.9%)
Size and weight: 132x140xh 60 mm - 450 g