

OLYMPIC 180 HOIST



Product Code:

HLA034

Description:

- The Olympic 180 Hoist benefits from the rigorous ergonomic design procedures and quality checks.
- It has been constructed from high quality powder coated steel that will provide years of reliable service.
- The hoist boasts an impressive 180kg lifting capacity, powered leg opening and market leading Linak® actuators and motors.

Safe Working Load:

185kg (29 Stone)

Lifting Range:

1755mm

Product Length:

1265mm

Product Height:

1335mm/2010mm (highest position/lowest position)

Leg Inside Width:

490mm/1075mm (closed/open)

Castor Size:

75mm/100mm (Front/Rear)

Carrier Bar Type

3 Hooked safety attachments for loop slings

Key Features:

- Cost effective solution without compromising on quality
- Robust powder coated steel frame
- Linak control system
- Electric leg opening as standard
- Red emergency lowering screw
- Easy to position with 4 large swivel castors and ergonomic handle
- Tool free assembly
- Single person installation
- Can lift a patient from the floor

Warranty:

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Hoist is guaranteed for five years from date of purchase.