



Pedal de elevação hidráulica dos 2 lados / Hydraulic two-sided lifting pedal
Pédale de levage hydraulique à deux côtés / Pedal de elevación hidráulica de los 2 lados

REF.	10.MA.1795
Descrição / Description Description / Descripción	Maca de Banho c/ sistema de elevação hidráulico Shower Trolley with hydraulic lifting system Chariot de douche avec système de levage hydraulique Carro de ducha con sistema de elevación hidráulico
Material Artigo / Material Article Materiél Article / Material artículo	Aço com tratamento de zincagem e acabamento a tinta poliéster Steel with zinc coating treatment and polyester paint finishing Acier avec traitement de revêtement au zinc et finition de peinture polyester Acero con tratamiento de cincado y acabado a tinta poliéster
Dimensão do leito / Dimension of bed Dimension du lit / Dimensión del lecho	1800x600x200mm
Dimensões Artigo / Dimensions Article / Article Dimensions Dimensiones Artículo	600x1878x410-830mm
Peso máximo paciente / Maximum patient weight / Poids maximal patient Peso máximo usuário	170kg
Peso Líquido / Net Weight Poids Net / Peso Neto	73kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	95kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,831m ³

