

Deluxe Wireless Floor Sensor Mat



Picture for illustrative purposes only

Product Code:

RG16

Description:

This floor sensor mat works will trigger an alarm when someone applies pressure by stepping onto the mat.

Features:

- Triggers alarm when pressure is released
- Durable and long-lasting
- Care instructions printed on pad
- Versatile and economical
- 1 year guarantee

Works in conjunction with:

- RG2 - Ramblegard Wireless Monitor
- RG9 - Transmitter

Technical Specifications:

Cable: for 433MHz Cordless transmitter

Material: PVC & PET

Dimensions (LxWxD): 122 x 61cm

Net Weight: 2kg

Operating temp.: 5°C to 40°C

Operating weight: 10 - 200 kg

Recommended Cleaning: 1:100 solution of household bleach (5.25-6.15% sodium hypo-chlorite)

**It stands out from its competitors,
it is effective and it's economical.**

MEDGUARD
Professional Healthcare Supplies

Medguard Professional Healthcare

Tel: +353 (0) 1 835 2411

Fax: +353 (0) 1 969 5009

Email: service@medguard.ie

Website: www.medguard.ie